

# Radioactive and Non-radioactive DNA Labeling by Nick-translation

This protocol is for the Radioactive and Non-radioactive DNA Labeling by Nick-translation.

## I. Radioactive DNA labeling by nick-translation

1. Mix the following components:

10x reaction buffer for DNA Polymerase I	2.5 $\mu$ L
Mixture of 3 dNTPs, 1 mM* (without the labeled dNTP)	1.25 $\mu$ L
[alpha-32P]-dNTP, ~ 110 TBq/mmol (3000 Ci/mmol)	1.85-3.7 MBq (50-100 $\mu$ Ci)
Thermo Scientific DNase I, RNase-free (Cat #EN0521, EN0523, EN0525) freshly diluted to 0.002 U/ $\mu$ L**	1 $\mu$ L
Thermo Scientific DNA Polymerase I (Cat #EP0041, EP0042)	0.5-1.5 $\mu$ L (5-15 U)
Template DNA	0.25 $\mu$ g
Water, nuclease-free	to 25 $\mu$ L
<b>Total volume</b>	<b>to 10 <math>\mu</math>L</b>

2. Immediately incubate at 15 °C for 15-60 minutes.
3. Terminate the reaction by adding 1  $\mu$ L of 0.5 M EDTA, pH 8.0.
4. Take an aliquot (1  $\mu$ L) to determine the efficiency of the label incorporation. A specific activity of DNA at least  $10^8$  cpm/ $\mu$ g DNA is expected.
5. If needed, the labeled DNA may be separated from the unincorporated radioactive precursors on Sephadex G-50 or Bio-Gel P-60 column or using spin column (Thermo Scientific GeneJET PCR Purification Kit, (Cat #K0701).

### Note:

\* To prepare a mixture of three non-labeled dNTPs (1 mM of each), mix 1  $\mu$ L aliquots of stock solutions of each dNTP (100 mM, Thermo Scientific Cat #R0181) with 97  $\mu$ L of Water, nuclease-free. These dNTP mixes can be stored at -20 °C for further use. \*\* The DNase I, RNase-free can be diluted with the 1x reaction buffer for DNA Polymerase I.

- The reaction volumes can be scaled up or down providing that the final concentrations of the components (DNA, dNTPs, labeled dNTP) are as indicated in the protocol.
- Radioactive DNA probes with higher specific activities can be prepared using two radioactively labeled dNTPs simultaneously. In this case, the composition of the unlabeled dNTP mix should be adjusted accordingly.

## II. Non-radioactive DNA labeling by nick-translation

The protocol above can be used for non-radioactive labeling by nick-translation using biotin-11-dUTP, fluorescein-12-dUTP, DIG-dUTP or aminoallyl-dUTP:

- Normal dTTP is substituted for labeled-dUTP at a molar ratio of 1:3-1:5,
- Reaction time is prolonged to 1-2 hours.

[www.thermoscientific.com/onebio](http://www.thermoscientific.com/onebio)

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

### Europe

Customer Service  
cs.molbio.eu@thermofisher.com

Technical Support  
ts.molbio.eu@thermofisher.com

Tel 00800 222 00 888  
Fax 00800 222 00 889

### United States

Customer Service  
cs.molbio@thermofisher.com

Technical Support  
ts.molbio@thermofisher.com

Tel 800 235 9880  
Fax 800 292 6088

### Canada

Customer Service  
cs.molbio@thermofisher.com

Technical Support  
ts.molbio@thermofisher.com

Tel 800 340 9026  
Fax 800 472 8322

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific