

# Analysis of Bromhexine Hydrochloride Using a Solid Core C18 Column

Joanne Jones, Thermo Fisher Scientific, Runcorn, Cheshire, UK

## Key Words

Accucore C18, bromhexine hydrochloride, solid core, Core Enhanced Technology

## Abstract

This application note demonstrates the use of the Thermo Scientific™ Accucore™ C18 HPLC column for the analysis bromhexine hydrochloride.

## Introduction

Bromhexine is a mucolytic agent that reduces mucus viscosity and aids the cilia (small hairs) in the respiratory tract in removing excess mucus. This application note describes the analysis of bromhexine using Thermo Scientific Accucore C18 HPLC columns.

Accucore HPLC columns use Core Enhanced Technology™ to facilitate fast and highly efficient separations. The 2.6 µm diameter particles are not totally porous, but instead have a solid core and a porous outer layer. The optimized phase bonding creates a series of high-coverage, robust phases. The carbon loading of the Accucore C18 column provides high retention of non-polar analytes via a predominantly hydrophobic interaction mechanism. The tightly controlled 2.6 µm diameter of Accucore particles results in much lower backpressures than typically seen with sub-2 µm materials.



## Experimental Details

Sample Handling	Part Number
Fisher Scientific™ HPLC grade water	W/0106/17
Fisher Scientific HPLC grade methanol	M/4062/17
8 mm Standard Opening Screw Thread Vial Convenience Kit, 2 mL Clear Vial with Patch, Black Polypropylene Closure with Red PTFE/White Silicone Septa	60180-600

## Sample Preparation

A primary standard of bromhexine hydrochloride was prepared in methanol at 1 mg/mL.

The working standard contained 2 µg/mL bromhexine hydrochloride in water / methanol (90:10 v/v)

Separation Conditions	Part Number
Instrumentation:	Thermo Scientific Dionex™ UltiMate™ 3000 RSLC HPLC System
Column:	Accucore C18 2.6 µm, 50 x 2.1 mm 17126-052130
Mobile phase:	0.1% formic acid in water / 0.1% formic acid in methanol (90:10, v/v)
Flow rate:	0.8 mL/min
Column temperature:	40 °C
Injection volume:	1 µL
Injection wash solvent:	Water / methanol (90:10, v/v)

## Results

Good retention and peak shape of bromhexine hydrochloride is demonstrated in Figure 1 and Table 1.




Figure 1: Chromatogram for bromhexine hydrochloride analyzed on an Accucore C18 2.6 µm, 50 x 2.1 mm column

	Bromhexine
Retention time (minutes)	1.44
%RSD on retention time	1.01
Asymmetry	1.12

Table 1: Results obtained from an Accucore C18 column, based upon 6 replicate injections

## Conclusion

Analysis of bromhexine can be achieved using an Accucore C18 column in under 2 minutes with excellent peak shape and good retention time reproducibility.

[thermoscientific.com/accucore](http://thermoscientific.com/accucore)

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**USA and Canada** +1 800 332 3331  
**France** +33 (0)1 60 92 48 34  
**Germany** +49 (0) 2423 9431 20 or 21  
**United Kingdom** +44 (0)1928 534110  
**Japan** +81 3 5826 1615

**China** +86 21 68654588 +86 10 84193588  
+86 20 83145199 800 810 5118  
**India** +91 22 6742 9494 +91 27 1766 2352  
**Australia** 1 300 735 292 (free call domestic)  
**New Zealand** 0800 933 966 (free call domestic)  
**All Other Enquiries** +44 (0) 1928 534 050

**Technical Support**  
**North America** +1 800 332 3331  
**Outside North America** +44 (0) 1928 534 440

**Thermo**  
SCIENTIFIC  
Part of Thermo Fisher Scientific