

## ► Physical Properties Table

Resin	Max. Use Temp. °C	HDT Temp. °C <sup>1</sup>	Brittleness Temp. °C <sup>13</sup>	Melting Point (°C)	Glass Transition Temp. (°C) <sup>20</sup>	Clarity	Microwavability <sup>2</sup>
<b>ABS</b>	86 to 90	75 to 98	0	210 to 230	105 to 109	Transparent to Opaque	--
<b>ACL</b>	120	177	-55	210 to 220	-83 & 100	Opaque	Marginal <sup>3</sup>
<b>COC</b>	70	70	-40	N/A	90	Transparent	Marginal <sup>3</sup>
<b>ECTFE</b>	150	116	<-76	242	85	Translucent	Yes
<b>EPDM</b>	146	<20	-68	N/A	-45	Opaque	Yes
<b>EPR</b>	145	< 20	-90	N/A	-54	Opaque	Yes
<b>ETFE</b>	150	104	-105	265	N/A	Translucent	Yes
<b>FEP</b>	205	70	-270	275	N/A	Translucent	Marginal <sup>3</sup>
<b>FEP 890</b>	205	70	-75	275	N/A	Transparent	Marginal <sup>3</sup>
<b>FLPE</b>	120	65	-100	125 to 138	-25	Translucent	No
<b>HDPE</b>	120	65	-100	125 to 138	-25	Translucent	No
<b>HIPS</b>	90	87	20	N/A	130	Opaque	Marginal <sup>3</sup>
<b>LDPE</b>	80	45	-100	85 to 125	-25	Translucent	Yes
<b>LLDPE 489 Tubing</b>	79	47	-73	85 to 125	-25	Translucent	Marginal <sup>3</sup>
<b>Modified PPE</b>	121	138	-40	240 to 320	164	Opaque	--
<b>PBS</b>	90	72	20	N/A	87-93	Transparent	Marginal <sup>3</sup>
<b>PC</b>	135	138	-135	N/A	154	Transparent	Marginal <sup>3</sup>
<b>PEI</b>	171	210	<20	N/A	215	Transparent Amber	Yes
<b>Permanox®</b>	174	127	-10	233	N/A	Transparent	Yes
<b>PES</b>	180	>200	-100	N/A	350 to 390	Transparent Amber	--
<b>PET</b>	65	21	-40	240 to 275	70	Transparent	Marginal <sup>3</sup>
<b>PETG</b>	70	70	-40	265	81	Transparent	Marginal <sup>3</sup>
<b>PFA</b>	260	73	-270	302 to 310	-80 & 90	Transparent	Yes
<b>PFA 870 Tubing</b>	260	73	-268	302 to 310	-80 & 90	Transparent	Yes
<b>PMMA</b>	50	79 - 107	20	85 to 105	N/A	Transparent	No
<b>PMP</b>	153 - 174	93 - 127	0 to 20	235	N/A	Transparent	Yes
<b>PP</b>	135	107	0	160 to 176	-20 to -5	Translucent	Yes
<b>PP 689 Tubing</b>	121	100	-4	160 to 176	-20 to -5	Translucent	Yes
<b>PPCO</b>	121	90	-40	150 to 175	-20	Translucent	Marginal <sup>3</sup>
<b>PS</b>	90	96	20	N/A	74 to 110	Transparent	No
<b>PSF</b>	165	174	-100	N/A	185 to 195	Transparent	Yes
<b>PTFE/TFE</b>	260	121	-100	320 to 330	120 to 130	Opaque	Yes
<b>PURS 280 Tubing</b>	85	32	-73	87 to 140	-30 to 0	Transparent Amber	No
<b>PVC 180 Tubing</b>	71	-32	-32	N/A	75 to 105	Transparent	Yes <sup>19</sup>
<b>PVC 380 Tubing</b>	79	-32	-32	N/A	75 to 105	Transparent	Yes <sup>19</sup>
<b>PVC 980 Tubing</b>	82	-32	-21	N/A	75 to 105	Transparent	No
<b>PVDF</b>	150	139	-62	141 to 178	-60 to -20	Translucent	Marginal <sup>3</sup>
<b>ResMer™</b>	130 to 150	200 to 300	20	200 to 270	90 to 110	Opaque	Marginal <sup>3</sup>
<b>SAN</b>	85	98	20	108	N/A	Transparent	Yes
<b>SILI 50/65 Tubing</b>	232	-46	-62	N/A	-130 to -120	Translucent	Yes
<b>Silicone Gaskets</b>	204	-46	-68	N/A	-130 to -120	Transparent to Opaque	Yes
<b>Thermanox®</b>	65	38	-60	220 to 240	80	Transparent	Marginal <sup>3</sup>
<b>TPE</b>	121	<23	<-50	N/A	N/A	Transparent to Opaque	Yes
<b>XLPE</b>	65	59	-118	N/A	N/A	Translucent	No
<b>Tritan™</b>	95	99	-40	282	108	Transparent	Marginal <sup>3</sup>

For abbreviations and footnotes, please see page 626-627.

## ► Physical Properties Table, continued

Resin	Sterilization <sup>4</sup>					Specific Gravity	Flexibility
	Autoclaving	EtO Gas	Dry Heat	Radiation	Disinfectants		
<b>ABS</b>	No	Yes	No	Yes	Some	1.05 to 1.16	Rigid
<b>ACL</b>	Yes <sup>5</sup>	Yes	No	No	Some	1.42 to 1.61	Rigid
<b>COC</b>	Yes	Yes	No	Yes	Some	1.03	Rigid
<b>ECTFE</b>	Yes	Yes	No	Yes	Yes	1.68	Rigid
<b>EPDM</b>	Yes	Yes	No	No	Some	0.86	Moderate to Excellent
<b>EPR</b>	Yes	Yes	No	No	Some	0.86	Moderate to Excellent
<b>ETFE</b>	Yes	Yes	No	Yes	Yes	1.7	Rigid
<b>FEP</b>	Yes	Yes	Yes	No	Yes	2.15	Excellent
<b>FEP 890</b>	Yes	Yes	Yes	No	Yes	2.17	Rigid
<b>FLPE</b>	No	Yes	No	Yes	Yes	0.95	Moderate
<b>HDPE</b>	No	Yes	No	Yes	Yes	0.95	Moderate
<b>HIPS</b>	No	Yes	No	Yes	Some	1.04	Rigid
<b>LDPE</b>	No	Yes	No	Yes	Yes	0.92	Good
<b>LLDPE 489 Tubing</b>	No	Yes	No	No	Some	0.92	Moderate
<b>Modified PPE</b>	Yes	Yes	No	Yes	Some	1.08	Rigid
<b>PBS</b>	No	Yes	No	Yes	Some	0.91 to 1.01	Rigid
<b>PC</b>	Yes <sup>5</sup>	Yes	No	Yes	Some	1.20	Rigid
<b>PEI</b>	Yes	Yes	Yes	Yes	Yes	1.28	Rigid
<b>Permanox®</b>	Yes	Yes	Yes	No	Yes	0.84	Rigid
<b>PES</b>	Yes	Yes	Yes	Yes	Some	1.37	Rigid
<b>PET</b>	No	Yes	No	Yes	Some	1.33	Moderate
<b>PETG</b>	No	Yes	No	Yes	Some	1.27	Moderate
<b>PFA</b>	Yes	Yes	Yes	No	Yes	2.17	Excellent
<b>PFA 870 Tubing</b>	Yes	Yes	Yes	No	Yes	2.15	Rigid
<b>PMMA</b>	No	No	No	Yes	Some	1.19	Rigid
<b>PMP</b>	Yes	Yes	Marginal	No	Yes	0.835	Rigid
<b>PP</b>	Yes	Yes	No	No	Yes	0.9	Rigid
<b>PP 689 Tubing</b>	Yes	Yes	No	No	Some	0.9	Rigid
<b>PPCO</b>	Yes	Yes	No	No	Yes	0.9	Moderate
<b>PS</b>	No	Yes	No	Yes	Some	1.05	Rigid
<b>PSF</b>	Yes	Yes	Yes	Yes	Yes	1.24	Rigid
<b>PTFE/TFE</b>	Yes	Yes	Yes	No	Yes	2.2	Rigid
<b>PURS 280 Tubing</b>	No	Yes	No	No	Some	1.18	Moderate
<b>PVC 180 Tubing</b>	Yes <sup>14</sup>	Yes	No	No	Some	1.19	Excellent
<b>PVC 380 Tubing</b>	Yes <sup>14</sup>	Yes	No	No	Some	1.2	Excellent
<b>PVC 980 Tubing</b>	No	Yes	No	No	Some	1.2	Moderate
<b>PVDF</b>	Yes	Yes	No	Yes	Yes	1.78	Rigid
<b>ResMer™</b>	Yes	Yes	No	Yes	Some	1.15 to 1.50	Rigid
<b>SAN</b>	No	Yes	No	Yes	No	1.07	Rigid
<b>SILI 50/65 Tubing</b>	Yes	Yes	Yes	Yes	Some	1.15	Excellent
<b>Silicone Gaskets</b>	Yes	Yes	Yes	Yes	Yes	1.1 to 1.5	Moderate to Excellent
<b>Thermanox®</b>	No	Yes	No	Yes	Some	1.3	Moderate
<b>TPE</b>	Yes	Yes	No	Yes	Some	0.9	Moderate to Excellent
<b>XLPE</b>	No	Yes	No	Yes	Yes	0.93	Rigid
<b>Tritan™</b>	No	Yes	No	Yes	Some	1.18	Moderate

For abbreviations and footnotes, please see page 626-627.

► **Physical Properties Table, continued**

Resin	Permeability (cc.-mil/100in <sup>2</sup> -24 hr.-atm)			Permeability (cc.-mm/m <sup>2</sup> -24 hr.-Bar)			Water Vapor Transmission Rate		Water Adsorption (%)
	N <sub>2</sub>	O <sub>2</sub>	CO <sub>2</sub>	N <sub>2</sub>	O <sub>2</sub>	CO <sub>2</sub>	(grams/m <sup>2</sup> ) <sup>15</sup>	(grams/100 in <sup>2</sup> ) <sup>16</sup>	
ABS	--	--	--	--	--	--	--	--	0.3 to 0.7
ACL	0.5	1.0	4.1	0.2	0.4	1.6	--	--	0.6
COC	15,871	5,032	1,742	6,167	1,955	677	40.00	2.58	<0.01
ECTFE	10	25	110	3.9	9.7	42.7	3.2	0.2	<0.1
EPDM	25 - 150	75 - 650	800 - 8000	9.7 - 58	29 - 253	3114 - 3108	--	--	0.05
EPR	25 - 150	75 - 650	800 - 8000	9.7 - 58	29 - 253	3114 - 3108	--	--	0.05
ETFE	30	100	250	12	39	97	1.65	0.11	0.03
FEP	320	750	2,200	124.34	291.41	854.82	6.200	0.40	<0.01
FEP 890	334	1,000	2,251	129.62	388.55	874.63	6.2	0.40	<0.01
FLPE	42	185	580	16	72	225	4.6	0.30	<0.01
HDPE	42	185	580	16	72	225	4.6 - 6.2	0.30 - 0.40	<0.01
HIPS	20 - 25	300 - 400	1,000 - 1,500	7.7 - 9.7	116 - 155	388 - 582	1220 - 6102	79 - 394	0.07
LDPE	180	500	2,700	70	194	1,049	15 - 23	1.0 - 1.5	<0.01
LLDPE 489 Tubing	334	1,001	4,670	130	389	1,815	15 - 23	1.0 - 1.5	<0.01
Modified PPE	--	1,000	--	--	389	--	69.70	4.50	0.10 to 0.15
PBS	20 - 25	300 - 400	1,000 - 1,500	7.8 - 9.7	116 - 155	388 - 582	2.20	0.14	0.07 to 0.09
PC	50	300	1,075	19.43	116.57	418	115.00	7.42	0.35
PEI	19	37	171	7.2	14.4	67	5.8	0.37	0.25
Permanox®	8,000	32,000	115,000	3,108	12,434	44,684	6.00	0.39	<0.01
PES	--	--	--	--	--	--	--	--	0.55
PET	0.7 - 1.0	13	70	0.27 - 0.39	5.1	27	18.13	1.17	0.25
PETG	10	25	125	3.9	9.7	49	18.13	1.17	0.13
PFA	291	881	2,260	113	342	878	2.0	0.13	<0.02
PFA 870 Tubing	300	1,084	2,502	117	421	972	2.0	0.13	<0.02
PMMA	2.78	12.4	68	1.08	4.8	26	55.20	3.56	0.35
PMP	8,000	32,000	115,000	3,108	12,434	44,684	775.00	50.00	0.01
PP	48	240	800	18.7	93	311	3.90	0.25	<0.02
PP 689 Tubing	67	417	1,501	26	162	583	3.9	0.25	<0.02
PPCO	45	200	650	17	78	253	4.40	0.28	<0.02
PS	20 - 25	300 - 400	1,000 - 1,500	7.8 - 9.7	116 - 155	388 - 582	1220 - 6102	79 - 394	0.05
PSF	55	300	700	21.40	117	272	--	--	0.3
PTFE/TFE	--	308	--	--	120	--	4.00	0.26	<0.01
PURS 280 Tubing	5.0 - 83	16.7 - 167	66.7 - 417	1.9 - 32.4	6.5 - 65	26 - 162	--	--	1.12
PVC 180 Tubing	8.3 - 33	16.7 - 100	167 - 584	3.2 - 13	2.6 - 39	64.8 - 227	15 - 80	1.0 - 5.2	0.15 to 0.75
PVC 380 Tubing	8.3 - 33	16.7 - 100	167 - 584	3.2 - 13	2.6 - 39	64.8 - 227	15 - 80	1.0 - 5.2	0.15 to 0.75
PVC 980 Tubing	8.3 - 33	16.7 - 100	167 - 584	3.2 - 13	2.6 - 39	64.8 - 227	15 - 80	1.0 - 5.2	0.15 to 0.75
PVDF	9	14	505	3.5	5.4	196	29.76	1.9	0.05
ResMer™	--	--	--	--	--	--	--	--	0.01
SAN	--	--	--	--	--	--	--	--	0.2
SILI 50/65 Tubing	46,116	132,762	335,741	17,918	51,585	130,453	--	--	0.1
Silicone Gaskets	--	--	--	--	--	--	--	--	0.1
Thermanox®	1.60	8.1	38.7	0.62	3.15	15.04	--	--	0.25
TPE	31 - 145	85 - 646	900 - 8634	12.05 - 56	33.03 - 251	350 - 3355	--	--	0.05 to 0.1
XLPE	42	185	580	16	72	225	4.6 - 6.2	0.3 - 0.4	<0.01
Tritan™	--	--	--	--	--	--	--	--	0.43

For abbreviations and footnotes, please see page 626-627.

## ► Physical Properties Table, continued

Resin	Non-Cytotoxicity <sup>6</sup>	Suitability for Food and Bev. Use <sup>7</sup>	Regulation Part 21 CFR	Refractive Index	Transparency (%)	Haze (%)	Hardness (Rockwell R/M or Shore D/A)	Coefficient of Friction (static)
<b>ABS</b>	--	--	--	--	78 to 90	2.5 to 11	93 to 115 R	--
<b>ACL</b>	Yes	Yes <sup>8</sup>	177.2480	1.48 to 1.51	0	--	94 - 120 R	0.07-0.20
<b>COC</b>	Yes	Yes	177-1520	1.533	91	--	--	--
<b>ECTFE</b>	Yes	Yes	177.1380	1.4200	--	--	90 R / 75 D	0.19
<b>EPDM</b>	--	Yes <sup>9</sup>	177.2600	--	0	--	45 - 87 A	0.725
<b>EPR</b>	--	Yes <sup>9</sup>	177.2600	--	0	--	45 - 87 A	0.853
<b>ETFE</b>	Yes	Yes	177.1380	1.3580	--	--	50 R	--
<b>FEP</b>	Yes	Yes	177.1550	1.341 to 1.347	96	--	60 D / 25 to 45 R	0.03 to 0.40
<b>FEP 890</b>	--	Yes	177.1550	1.347	96	--	58 D	--
<b>FLPE</b>	Yes	Yes <sup>9</sup>	177.1615	1.5400	--	--	60 to 73 D	0.29
<b>HDPE</b>	Yes	Yes <sup>9</sup>	177.1520	1.5400	--	--	60 to 73 D	0.29
<b>HIPS</b>	Yes	Yes	177.1640	--	0	--	--	--
<b>LDPE</b>	Yes	Yes <sup>9</sup>	177.1520	1.5100	--	--	41 to 46 D	--
<b>LLDPE 489 Tubing</b>	--	Yes <sup>9</sup>	177.1520	1.5100	--	--	50 D	--
<b>Modified PPE</b>	Yes	Yes	177.2460 <sup>17</sup>	--	0	--	115 R to 120 R / 78 M	--
<b>PBS</b>	Yes	Yes	177.1640	1.535 to 1.573	90 to 95	--	53 to 64 D	--
<b>PC</b>	Yes	Yes	177.1580	1.5860	91	--	70 M	0.31
<b>PEI</b>	--	Yes	177.1595	1.4600	58	--	125 R	--
<b>Permanox®</b>	Yes	Yes	177.1520	1.46	94	--	87 R	--
<b>PES</b>	Yes	Yes	177.2440	1.65	90	--	120 to 127 R	88 M
<b>PET</b>	Yes	Yes	177.1315	1.5750	85	--	112 R	0.2 to 0.4
<b>PETG</b>	Yes	Yes <sup>10</sup>	177.1315	1.57	91	--	106 R	--
<b>PFA</b>	Yes	No	--	1.3580	93.5	--	55 D	0.18
<b>PFA 870 Tubing</b>	--	Yes	177.1550	1.3580	93.5	--	60 D	--
<b>PMMA</b>	Yes	Yes	177.1010	1.48 to 1.57	92	--	92 to 100 M	--
<b>PMP</b>	Yes	Yes <sup>11</sup>	177.1520	1.4630	93	2.1	85 R	--
<b>PP</b>	Yes	Yes	177.1520	1.4735 to 1.490	--	--	80-110 R / 70-75 D	0.24
<b>PP 689 Tubing</b>	--	Yes	177.1520	1.4900	--	--	75 D	--
<b>PPCO</b>	Yes	Yes	177.1520	1.4735 to 1.5100	--	--	80-100 R / 70-75 D	0.24
<b>PS</b>	Yes	Yes	177.1640	1.5894	90	--	60 to 90 M	--
<b>PSF</b>	Yes	Yes	177.1655	1.6330	86	--	120 R	0.29
<b>PTFE/TFE</b>	Yes	Yes	177.1550	1.3500	0	--	58 R / 50-65 D	0.04
<b>PURS 280 Tubing</b>	--	No	--	1.5 to 1.6	--	--	85 A	--
<b>PVC 180 Tubing</b>	--	Yes <sup>12</sup>	176.180 / 175.300	1.53 to 1.55	--	--	55 A	--
<b>PVC 380 Tubing</b>	--	Yes <sup>12</sup>	176.180 / 175.300	1.53 to 1.55	--	--	65 A	--
<b>PVC 980 Tubing</b>	--	Yes <sup>12</sup>	176.180 / 175.300	1.53 to 1.55	--	--	65 A	--
<b>PVDF</b>	Yes	Yes	177.2510	1.4200	--	--	78 D	0.3
<b>ResMer™</b>	--	--	--	--	--	--	--	--
<b>SAN</b>	--	Yes	177.1040	1.5700	89	--	125 R	--
<b>SILI 50/65 Tubing</b>	Yes	Yes <sup>9</sup>	177.2600	1.4035	--	--	50 A / 65 A	--
<b>Silicone Gaskets</b>	--	Yes <sup>9</sup>	177.2600	--	0 to 95	--	--	--
<b>Thermanox®</b>	Yes	Yes	177.1315	1.57	91	--	94 M	0.2 to 0.4
<b>TPE</b>	Yes	Yes	177.2600	--	0 to 95	--	--	--
<b>XLPE</b>	Yes	No	--	1.5400	--	--	50 to 70 R	0.1 to 0.2
<b>Tritan™</b>	Yes	Yes	174.5(d)(5) FCN 729	--	90	1	112 R	--

For abbreviations and footnotes, please see page 626-627.

## ► Physical Properties Abbreviations

**Table of Abbreviations**

Abbreviation	Name
ACL	Acetal
COC	Cyclic Olefin Copolymer
ECTFE	Ethylene Chlorotrifluoroethylene
EPDM	Ethylene Propylene Diene Monomer Rubber
EPR	Ethylene Propylene Rubber
ETFE	Ethylene Tetrafluoroethylene
FEP	Fluorinated Ethylene Propylene
FLPE	Fluorinated High-density Polyethylene
HDPE	High-density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
PC	Polycarbonate
PEI	Poetherimide
PET	Polyethylene terephthalate
PETG	Polyethylene terephthalate copolymer
PFA	Perfluoroalkoxy
PMMA	Polymethylmethacrylate
PMP	Polymethylpentene
PP	Polypropylene
PPCO	Polypropylene Copolymer
PPE	Modified Polyphenylene Ether
PS	Polystyrene
PSF	Polysulfone
PURS	Thermoplastic Polyester Polyurethane
PVDF	Polyvinylidenedifluoride
ResMer™	ResMer™ - Glass Filled
SAN	(No Suggestions)
Silicone Gaskets	Silicone - Gaskets, filled - typically red
Silicone 50/65 Tubing	Silicone - Tubing, and unfilled (translucent) o-rings
PTFE / TFE	Tetrafluoroethylene
TPE	Thermoplastic Elastomer
Thermanox®	Proprietary NUNC Resin
Permanox®	Proprietary NUNC Resin
Tritan™	Proprietary CoPolyester

## ► Physical Properties Footnotes

- 1 Heat Deflection Temperature is the temperature at which an injection molded bar deflects 0.01 in. when placed under 66 psig (ASTM D648) of pressure. Materials may be used above Heat Deflection Temperatures in non-stress applications; see Max. Use Temp.
- 2 Ratings based on five minute tests using 600 watts of power on exposed, empty labware. CAUTION: Do not exceed Max. Use Temp., or expose labware to chemicals which heating cause to attack the plastic or be rapidly absorbed.
- 3 The plastic will absorb and retain significant amounts of heat resulting in an unexpectedly hot surface.
- 4 **STERILIZATION:**  
Autoclaving (121°C, 15 psig for 20 minutes) – Clean and rinse items with distilled before autoclaving. (Always completely disengage thread before autoclaving.) Certain chemicals which have no appreciable effect on resins at room temperature may cause deterioration at autoclaving temperatures unless removed with distilled water before hand.  
EtO Gas – Ethylene Oxide: 100% EtO, EtO:Nitrogen mixture, EtO:HCFC mixture  
Dry Heat – exposure to 160°C for 120 minutes without stress/load on the polymer parts  
Disinfectants – Benzalkonium chloride, formalin/formaldehyde, hydrogen peroxide, ethanol, etc.  
Radiation – gamma or beta irradiation at 25 kGy (2.5 MRad) with unstabilized plastic.
- 5 Sterilizing reduces mechanical strength. Do not use PC vessels for vacuum applications if they have been autoclaved. Refer to Use and Care Guidelines for Nalgene™ Labware, for detailed information on sterilizing.
- 6 “Yes” indicates the resin has been determined to be non-cytotoxic, based on USP and ASTM biocompatibility testing standards utilizing an MEM elution technique with WI38 human diploid lung cell line.
- 7 Resins meet requirements of CFR21 section of Food Additives Amendment of the Federal Food and Drug Act. End users are responsible for validation of compliance for specific containers used in conjunction with their particular applications.
- 8 Acceptable for aqueous foods only, at temperatures up to 121°C / 250°F. Not sanctioned for use with alcoholic or fatty foods at any temperature.
- 9 Acceptable for:  
Nonacid, aqueous products; may contain salt, sugar or both (pH above 5.0)  
Dairy products and modifications; oil-in-water emulsions, high or low fat  
Moist bakery products with surface containing no free fat or oil  
Dry solids with the surfaces containing no free fat or oil (no end-test required) and under all conditions as described in Table 2 of FDA Regulation 177.1520 except condition A - high temperature sterilization (e.g. over 100°C / 212°F)
- 10 Acceptable for:  
Alcoholic foods containing not more than 15% (by volume) alcohol; fill and storage temperature not to exceed 49°C (120°F)  
Non-alcoholic foods of hot fill to not exceed 82°C (180°F) and 49°C (120°F) in storage.  
Not suitable for carbonated beverages or beer or packaging food requiring thermal processing.
- 11 Straight-sided jars, beakers and graduated cylinders only.
- 12 Acceptable for aqueous, oil, dairy, acidic, and alcoholic foods up to 71°C/160°F.
- 13 The brittleness temperature is the temperature at which an item made from the resin may break or cracked if dropped. This is not the lowest use temperature if care is exercised in use and handling.
- 14 The tubing will become opaque from absorbed water, see the current NALGENE® Labware catalog information on Autoclaving PVC Tubing for details.
- 15 WVTR = Water Vapor Transmission rate in g-mm/m<sup>2</sup> - 24 hr. - 1 Bar at 37°C and 90% Relative Humidity.
- 16 WVTR = Water Vapor Transmission rate in g-mm/100 in<sup>2</sup> - 24 hr. - 1 Bar at 37°C and 90% Relative Humidity.
- 17 Acceptable for:  
Frozen or refrigerated storage: Ready-prepared foods intended to be reheated in container at time of use  
1. Aqueous or oil-in-water emulsion of high- or low-fat.  
2. Aqueous, high- or low-free oil or fat.
- 18 Will vary based on exact composition
- 19 If microwaved in the presence of water; the tubing will become opaque from absorbed moisture, see the current Nalgene Labware catalog information on Autoclaving PVC Tubing for details.
- 20 Glass Transition Temperature is the reversible transition in amorphous or semi-amorphous materials from a hard and relatively brittle state into a flexible or rubber-like state without becoming a true liquid.