

Quantitative Analysis of Immunosuppressant Drugs in Whole Blood Using High Throughput LC-MS/MS

Marta Kozak¹, Raimund Ruzicka², Vani Bodepudi², Jeff Zonderman³, ¹Thermo Fisher Scientific, San Jose, CA, USA, ²Thermo Fisher Scientific, Fremont, CA, USA, ³Thermo Fisher Scientific, Franklin, MA, USA

Key Words

- TSQ Quantum Ultra
- Transcend LX-2 System
- Clinical Research
- Multiplexing

Introduction

Immunosuppressant drugs inhibit the body's immune system and are used in organ transplant patients to prevent organ rejection. Liquid chromatography-mass spectrometry (LC-MS/MS) is a widely-accepted technique for the determination of immunosuppressant drugs in whole blood by research laboratories.

This application note describes a fast, sensitive, reliable, and accurate LC-MS/MS quantitative method for use by research laboratories for the simultaneous analysis of tacrolimus, sirolimus, everolimus and cyclosporin A in whole blood.

Experimental Conditions

Sample Preparation

A protein precipitation solution was prepared by mixing MeOH containing internal standards (Ascomycin and Cyclosporin D) with ZnSO₄ solution. Blood samples were processed by adding precipitation solution. The mixture was vortexed and centrifuged. Supernatant was injected into the LC-MS/MS system.

HPLC

HPLC analysis was performed using Thermo Scientific Transcend LX-2 advanced multiplexing system. Samples were injected into a Thermo Scientific Javelin C18 guard column at 80 °C and analyzed with a 2-minute gradient method.

Mass Spectrometry

MS analysis was performed using a Thermo Scientific TSQ Quantum Ultra triple stage quadrupole mass spectrometer equipped with an atmospheric pressure chemical ionization (APCI) source in selective reaction monitoring (SRM) data acquisition mode. Optimized SRM parameters for all of the analytes and internal standards are shown in Table 1.

Table 1: Optimized SRM parameters

Compound	Parent Ion	Fragment Ion	Collision Energy	Tube Lens Offset
Tacrolimus	821.4	768.3	18	190
Sirolimus	931.6	864.5	15	190
Everolimus	975.7	908.4	16	190
Ascomycin	809.4	756.4	18	190
Cyclosporin A	1219.9	1202.9	17	190
Cyclosporin D	1234.0	1216.9	17	190

Results and Discussion

Figure 1 displays the representative limits of quantitation (LOQ) chromatograms for tacrolimus, sirolimus, everolimus, cyclosporin A, and the internal standards. As shown in Tables 2 and 3, the intra- and inter-day variability were




Figure 1:
Chromatogram of
lowest calibration
standard

excellent. For each analyte, intra-day variability was determined by processing and analyzing 5 replicates of each QC sample. Inter-day variability was determined with 5 replicates of each QC sample in 3 different batches. The method tested negatively for all interferences and cross-reactivity. No ion suppression or enhancement was observed.

Assay performance summary

Target Analytes	Tacrolimus, Sirolimus, Everolimus and Cyclosporin A
Matrix	Whole blood
LOQ	1 ng/mL (Tacrolimus, Sirolimus, Everolimus) 10 ng/mL (Cyclosporin A)
Assay Linearity	1-50 ng/mL (Tacrolimus, Sirolimus, Everolimus) 10-2000 ng/mL (Cyclosporin A)
Analysis Time	2.0 min; 1.0 min with column multiplexing

Table 2: Intra-day variability (%RSD)

Analyte	QC1	QC2	QC3	QC4	QC5
Tacrolimus	6.8	4.6	4.9	-	-
Sirolimus	6.7	6.1	3.9	-	-
Everolimus	8.6	5.1	4.5	-	-
Cyclosporin A	5.5	4.6	3.4	3.4	3.5

Table 3: Inter-day variability (%RSD)

Analyte	QC1	QC2	QC3	QC4	QC5
Tacrolimus	4.2	4.1	1.7	-	-
Sirolimus	4.4	7.0	7.5	-	-
Everolimus	7.5	2.3	6.8	-	-
Cyclosporin A	1.8	2.0	2.4	1.7	4.7

Conclusion

A fast, sensitive, reliable and accurate method was developed for the quantification of tacrolimus, sirolimus, everolimus and cyclosporin A in whole blood by research laboratories. The use of column multiplexing technology allows for a 1 min analytical method, which enhances sample throughput.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Africa-Other

+27 11 570 1840

Australia

+61 2 8844 9500

Austria

+43 1 333 50 34 0

Belgium

+32 2 482 30 30

Canada

+1 800 530 8447

China

+86 10 8419 3588

Denmark

+45 70 23 62 60

Europe-Other

+43 1 333 50 34 0

**Finland/Norway/
Sweden**

+46 8 556 468 00

France

+33 1 60 92 48 00

Germany

+49 6103 408 1014

India

+91 22 6742 9434

Italy

+39 02 950 591

Japan

+81 45 453 9100

Latin America

+1 608 276 5659

Middle East

+43 1 333 50 34 0

Netherlands

+31 76 579 55 55

South Africa

+27 11 570 1840

Spain

+34 914 845 965

Switzerland

+41 61 716 77 00

UK

+44 1442 233555

USA

+1 800 532 4752

www.thermo.com

Legal Notices

©2010 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific Inc. products. It is not intended to encourage use of these products in any manners that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

View additional Thermo Scientific LC/MS application notes at: www.thermo.com/appnotes



Thermo Fisher Scientific,
San Jose, CA USA is ISO Certified.

AN62944_E 12/09M