

Determination of Ultratrace Elements in Semiconductor Grade Hydrochloric Acid using the Thermo Scientific iCAP Qs ICP-MS

¹Tomoko Vincent, ²Tokutaka Ikemoto, ²Tomoko Aoki, ²Makiko Shimura

¹Thermo Fisher Scientific, Bremen, Germany, ²Thermo Fisher Scientific, Yokohama, Japan

Key Words

Hydrochloric Acid, Semiconductor, Cold Plasma, KED, CCT

Goal

To develop a method for the ultratrace determination of metals in semiconductor grade hydrochloric acid (HCl) using the Thermo Scientific iCAP Qs ICP-MS.

Introduction

The organization SEMI (Semiconductor Equipment and Materials International, www.semi.org) has announced a guideline (defined as “Future Chemical – Purity Needs”) for maximum trace metal concentrations to be as low as 10 ng·L⁻¹ in process chemicals. For example SEMI specification C27-0708 defines maximum metal contamination levels in concentrated hydrochloric acid (HCl).

HCl is used in silicon wafer cleaning processes to remove dust, organic contamination and any thin oxide layer. This cleaning process, called SC-2 (Standard Clean, 2nd step) is performed with a 1-15% HCl solution.

Due to a series of complex Cl and Ar based polyatomic interferences, ³⁹K, ⁵¹V, ⁵²Cr, ⁷⁴Ge and ⁷⁵As are often difficult to analyze at ultratrace concentrations in HCl.

Sample preparation

Pre-cleaned PFA bottles were used for the preparation of all blanks, standards and samples. Standards at concentrations of 10, 25, 50 and 100 ng·L⁻¹ were prepared gravimetrically by adding the appropriate quantity of a multi-elemental SPEX CertiPrep™ stock solution directly to concentrated HCl (32-35% HCl Optima™ grade from Fisher Chemical), and diluted with ultrapure water to make 20% HCl. 20% HCl solutions for analysis was used for the rinse and blank solutions. Spike tests were performed at 10 ng·L⁻¹.

Method

The instrument configuration and operation parameters are shown in Table 1. Please note that the Thermo Scientific™ iCAP™ Qs ICP-MS was not installed in a clean room.

Table 1: Instrument configuration and operation parameters.






Parameter	Value
Spray Chamber	Quartz cyclonic
Nebulizer	Microflow PFA-100 (self-aspirating)
Injector	2.0 mm I.D., sapphire
Interface	Cold plasma platinum sampler and skimmer
Extraction Lens	Cold plasma lens kit

Mode	RF Power	QCell Technique
CP	580 W	No Gas
CP-NH ₃	580 W	1% NH ₃ in 99% He, 7.0 mL·min ⁻¹
CCT-NH ₃	1550 W	50% NH ₃ in 50% He, 0.8 mL·min ⁻¹
CCT-O ₂	1550 W	50% O ₂ in 50% He, 2.0 mL·min ⁻¹
KED-He	1550 W	100% He, 5.0 mL·min ⁻¹
KED-NH ₃	1550 W	1% NH ₃ in 99% He, 4.0 mL·min ⁻¹

Results

As can be seen in Figures 1 to 5, the powerful collision cell capabilities of the QCell in the iCAP Qs ICP-MS provide the advanced performance required for the sensitive and accurate determination of difficult elements, such as K, V, Cr, Ge, As in 20% HCl.

- Cold plasma inhibits the formation of Cl and Ar based interferences (on K and Cr).
- Cold plasma promotes the formation of AsO⁺ (shifting the target analyte to an interference free mass range).
- The combination of automatic low mass cut off with different QCell modes (CCT, KED) effectively suppress Cl based interferences (on V and Ge).

Figure 1: Calibration curve for ^{39}K in 20% HCl (interference $^1\text{H}_2\ ^{37}\text{Cl}^+$).Figure 2: Calibration curve for ^{51}V in 20% HCl (interference $^{35}\text{Cl}^{16}\text{O}^+$).Figure 3: Calibration curve for ^{52}Cr in 20% HCl (interference $^{35}\text{Cl}^{16}\text{O}^+\text{H}^+$).Figure 4: Calibration curve for ^{74}Ge in 20% HCl (interference $^{37}\text{Cl}_2^+$).Figure 5: Calibration curve for ^{75}As (as $^{75}\text{As}^{16}\text{O}^+$) in 20% HCl (interference $^{40}\text{Ar}^{35}\text{Cl}^+$).

The capability of the iCAP Qs ICP-MS to reliably analyze 20% HCl is demonstrated by the excellent LoD and BEC value shown in Table 2. Spike recoveries (for only 10 $\text{ng}\cdot\text{L}^{-1}$) of between 91% to 115% for 29 elements in 20% HCl demonstrate the excellent interference removal and accuracy achieved by the iCAP Qs ICP-MS using the method described.

Table 2: iCAP Qs ICP-MS performance for 20% HCl. As ($^{75}\text{As}^{16}\text{O}^+$) A detection limit of $<0.001\ \text{ng}\cdot\text{L}^{-1}$ was recorded for ^{75}As since all repeats at $^{75}\text{As}^{16}\text{O}^+$ in the blank gave the same intensities "0".

Analyte	Mode	LoD ($\text{ng}\cdot\text{L}^{-1}$)	BEC ($\text{ng}\cdot\text{L}^{-1}$)	Recovery (%)
^7Li	CP	0.07	0.05	108
^{23}Na	CP	0.2	0.4	105
^{24}Mg	CP	0.3	0.2	110
^{27}Al	CP	0.7	0.8	115
^{39}K	CP-NH ₃	0.4	1.8	102
^{40}Ca	CP	1.4	9.5	104
^{51}V	CCT-NH ₃	0.9	3.2	93
^{52}Cr	CP-NH ₃	1.4	7.8	111
^{55}Mn	CP-NH ₃	0.5	0.4	105
^{57}Fe	CP	1.0	3.1	100
^{59}Co	CP	0.5	0.3	94
^{60}Ni	CP	2.6	4.7	95
^{63}Cu	CP	1.8	3.1	104
^{66}Zn	KED-NH ₃	3.0	9.6	96
^{71}Ga	KED-He	1.3	0.8	96
^{74}Ge	KED-NH ₃	3.4	4.1	99
^{75}As (AsO ⁺)	CP	<0.001	<0.001	99
^{80}Se	CCT- O ₂	0.3	15	93
^{85}Rb	CP	0.2	0.06	108
^{88}Sr	KED-He	0.2	0.02	105
^{111}Cd	KED-He	0.6	0.1	92
^{115}In	KED-He	0.2	0.3	93
^{133}Cs	KED-He	0.01	0.01	93
^{138}Ba	KED-He	0.01	0.01	98
^{202}Hg	KED-He	0.6	1.3	91
^{205}Tl	KED-He	0.2	0.3	92
^{208}Pb	KED-He	0.2	1.2	94
^{209}Bi	KED-He	0.4	1.6	93
^{238}U	KED-He	0.03	0.03	93

Conclusion

The Thermo Scientific iCAP Qs ICP-MS has been shown to offer the high sensitivity and freedom from interferences required for the measurement of ultratrace ($\text{ng}\cdot\text{L}^{-1}$) concentration levels in semiconductor grade HCl. Flexible operation combined with the high performance of the QCell provide significant improvements ideally suited for the requirements of advanced semiconductor applications.

www.thermofisher.com

©2016 Thermo Fisher Scientific Inc. All rights reserved. CertiPrep is a trademark of SPEX CERTIPREP, INC. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 9 3291 0200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00
Russia/CIS +43 1 333 50 34 0

Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand