

# Determination of Ultratrace Elements in Semiconductor Grade Sulfuric Acid using the Thermo Scientific iCAP Qs ICP-MS

Tomoko Vincent, Christoph Wehe  
Thermo Fisher Scientific, Bremen, Germany

## Key Words

Sulfuric Acid, Semiconductor, Cold Plasma, KED, CCT

## Goal

To develop a method for the ultratrace determination of metals in semiconductor grade sulfuric acid using the Thermo Scientific iCAP Qs ICP-MS.

## Introduction

Concentrated sulfuric acid ( $\text{H}_2\text{SO}_4$ ) is used in the semiconductor industry to remove organic substances from the surface of silicon wafers. Since  $\text{H}_2\text{SO}_4$  comes into contact with wafer surfaces, it is necessary to monitor its trace elemental impurities.

High sulfur matrices are problematic in ICP-MS analysis due to the formation of spectral interferences that are preferentially created due to their low ionization energies in the ICP ion source. The most challenging elements in the trace determination of sulfuric acid are Ti, V, Cr, Zn and Ge. All major isotopes of these metals are severely interfered by polyatomic species generated by the sample matrix (see Table 1).

Table 1: Typical target isotopes with commonly observed interferences.

Isotope	Abundance	Matrix-Based Interferences
$^{47}\text{Ti}$	7.3%	$^{33}\text{S}^{14}\text{N}^+$ , $^{32}\text{S}^{15}\text{N}^+$ , $^{32}\text{S}^{14}\text{N}^+\text{H}^+$
$^{48}\text{Ti}$	73.8%	$^{32}\text{S}^{16}\text{O}^+$
$^{51}\text{V}$	100%	$^{33}\text{S}^{18}\text{O}^+$ , $^{34}\text{S}^{16}\text{O}^+\text{H}^+$ , $^{32}\text{S}^{18}\text{O}^+\text{H}^+$
$^{52}\text{Cr}$	83.8%	$^{34}\text{S}^{18}\text{O}^+$ , $^{33}\text{S}^{18}\text{O}^+\text{H}^+$
$^{53}\text{Cr}$	9.5%	$^{34}\text{S}^{18}\text{O}^+\text{H}^+$
$^{64}\text{Zn}$	48.6%	$^{32}\text{S}_2^+$ , $^{32}\text{S}^{16}\text{O}_2^+$
$^{66}\text{Zn}$	27.9%	$^{32}\text{S}^{34}\text{S}^+$ , $^{33}\text{S}_2^+$ , $^{34}\text{S}^{16}\text{O}_2^+$ , $^{36}\text{S}^{16}\text{O}^{14}\text{N}^+$
$^{68}\text{Zn}$	18.8%	$^{34}\text{S}_2^+$ , $^{32}\text{S}^{36}\text{S}^+$ , $^{34}\text{S}^{16}\text{O}^{18}\text{O}^+$
$^{72}\text{Ge}$	27.7%	$^{40}\text{Ar}^{32}\text{S}^+$
$^{74}\text{Ge}$	35.5%	$^{40}\text{Ar}^{34}\text{S}^+$

## Sample Preparation

Pre-cleaned PFA bottles were used for the preparation of all blanks, standards and samples. Concentrated  $\text{H}_2\text{SO}_4$  (98% Optima grade  $\text{H}_2\text{SO}_4$ , from Fisher Chemical) was 10-fold diluted with ultrapure water before analysis. Standards at concentrations of 10, 25, 50 and 100  $\text{ng}\cdot\text{L}^{-1}$  were prepared gravimetrically by adding the appropriate quantity of a multi-elemental SPEX CertiPrep™ stock. 9.8%  $\text{H}_2\text{SO}_4$  was used for the rinse and blank solutions. Spike tests were performed at 10  $\text{ng}\cdot\text{L}^{-1}$ .

## Method

The instrument configuration and operation parameters are shown in Table 2. Please note, the Thermo Scientific™ iCAP™ Qs ICP-MS was not installed in a clean room.

Table 2: Instrument configuration and operation parameters.

Parameter	Value
Spray Chamber	Quartz cyclonic
Nebulizer	MicroFlow PFA-100 (self-aspirating)
Injector	2.0 mm I.D., Sapphire
Interface	Cold plasma platinum sampler and skimmer
Extraction Lens	Cold plasma lens kit

Mode	RF Power	QCell Technique
CP	580 W	-
CP-NH <sub>3</sub>	580 W	1% NH <sub>3</sub> in He, 7.0 mL·min <sup>-1</sup>
CCT-Cluster	1550 W	50% NH <sub>3</sub> in He, 0.8 mL·min <sup>-1</sup>
CCT-NH <sub>3</sub>	1550 W	1% NH <sub>3</sub> in He, 10.0 mL·min <sup>-1</sup>
KED-He	1550 W	100% He, 5.0 mL·min <sup>-1</sup>
KED-NH <sub>3</sub>	1550 W	1% NH <sub>3</sub> in He, 4.0 mL·min <sup>-1</sup>

## Results

Table 3 shows the performance of the analysis of 9.8%  $\text{H}_2\text{SO}_4$  with the iCAP Qs ICP-MS:

- Excellent semiconductor level LoD and BEC were produced for all analytes through a combination of different measurement modes.
- Outstanding performance was achieved for V, Cr, Zn and Ge indicating the effective suppression of all matrix induced interferences (see Figure 1 for V).
- CCT-Cluster mode (50%  $\text{NH}_3$  in He) effectively shifts the target analyte to an interference free region enabling ultratrace Ti quantification (see Figure 2).
- Accurate spike recoveries from 92% to 108% were obtained for 27 elements at  $10 \text{ ng}\cdot\text{L}^{-1}$ .
- Cold plasma is suitable for the analysis of mineral acid samples such as 9.8% sulfuric acid.




Figure 1: Calibration curve for  $^{51}\text{V}$  in 9.8%  $\text{H}_2\text{SO}_4$ .




Figure 2: Calibration curve for  $^{48}\text{Ti}^{14}\text{N}^1\text{H}^{(14}\text{N}^1\text{H}_3)_3^+$  in 9.8%  $\text{H}_2\text{SO}_4$ .

## Conclusion

The Thermo Scientific iCAP Qs ICP-MS provides excellent performance for the ultratrace determination of metals in semiconductor grade  $\text{H}_2\text{SO}_4$ . The flexible combination of different analysis modes has been shown to be ideally suited for the ultratrace metal determination in advanced semiconductor applications.

[www.thermofisher.com](http://www.thermofisher.com)

©2016 Thermo Fisher Scientific Inc. All rights reserved. CertiPrep is a trademark of SPEX CERTIPREP, INC. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Africa** +43 1 333 50 34 0  
**Australia** +61 3 9757 4300  
**Austria** +43 810 282 206  
**Belgium** +32 53 73 42 41  
**Canada** +1 800 530 8447  
**China** 800 810 5118 (free call domestic)  
 400 650 5118

**Denmark** +45 70 23 62 60  
**Europe-Other** +43 1 333 50 34 0  
**Finland** +358 9 3291 0200  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9494  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Latin America** +1 561 688 8700  
**Middle East** +43 1 333 50 34 0  
**Netherlands** +31 76 579 55 55  
**New Zealand** +64 9 980 6700  
**Norway** +46 8 556 468 00  
**Russia/CIS** +43 1 333 50 34 0

**Singapore** +65 6289 1190  
**Spain** +34 914 845 965  
**Sweden** +46 8 556 468 00  
**Switzerland** +41 61 716 77 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752

Table 3: iCAP Qs ICP-MS performance for the analysis of 9.8%  $\text{H}_2\text{SO}_4$ .

Analyte	Mode	LoD (ng·L <sup>-1</sup> )	BEC (ng·L <sup>-1</sup> )	Recovery (%)
<sup>7</sup> Li	CP	0.3	0.9	99
<sup>23</sup> Na	CP	0.7	2.6	107
<sup>24</sup> Mg	CP-NH <sub>3</sub>	0.9	0.8	102
<sup>27</sup> Al	CP	0.9	1.1	99
<sup>39</sup> K	CP	2.2	8.5	102
<sup>40</sup> Ca	CP	5.9	17.0	103
<sup>48</sup> TiNH(NH <sub>3</sub> ) <sub>3</sub>	CCT-Cluster	2.7	21.7	97
<sup>51</sup> V	CCT-NH <sub>3</sub>	0.6	0.2	100
<sup>53</sup> Cr	CP-NH <sub>3</sub>	2.7	1.7	103
<sup>55</sup> Mn	CP	0.3	0.2	102
<sup>56</sup> Fe	CP	1.0	1.8	100
<sup>58</sup> Ni	CP	1.5	2.1	98
<sup>59</sup> Co	CP	0.3	0.1	99
<sup>65</sup> Cu	CP	1.9	1.1	103
<sup>68</sup> Zn	KED-NH <sub>3</sub>	1.9	2.9	99
<sup>71</sup> Ga	CP	0.5	0.1	101
<sup>74</sup> Ge	KED-He	0.1	0.1	108
<sup>75</sup> As	KED-NH <sub>3</sub>	0.7	0.8	100
<sup>85</sup> Rb	KED-He	1.0	0.6	105
<sup>88</sup> Sr	KED-NH <sub>3</sub>	0.5	0.4	107
<sup>111</sup> Cd	KED-He	1.2	0.5	107
<sup>115</sup> In	KED-He	0.3	0.1	96
<sup>138</sup> Ba	KED-NH <sub>3</sub>	1.1	0.5	92
<sup>205</sup> Tl	KED-He	0.6	1.4	97
<sup>208</sup> Pb	KED-He	1.5	3.1	98
<sup>209</sup> Bi	KED-He	0.5	0.4	101
<sup>238</sup> U	KED-He	0.4	0.3	105