

Targeted Quantitative Protein Analysis in Human Serum, Using High Resolution Multiple Selected Reaction Monitoring Assays

Reiko Kiyonami, Scott Peterman, and Ken Miller, Thermo Fisher Scientific, San Jose, CA, USA

Key Words

- TSQ Quantum Ultra™
- Surveyor HPLC™
- Biomarkers
- H-mSRM
- Human serum complex sample
- Proteomics
- Targeted protein quantitation

Introduction

A common endpoint for a biomarker discovery experiment is a list of putative marker proteins. A reasonable next step is to then perform targeted quantitative measurements of these proteins in an expanded patient population to assess their validity as markers. Analytical accuracy and precision are required for unambiguous quantitative analysis of targeted proteins from very complex mixtures. Wide dynamic range and high sensitivity are critical for detecting low abundance proteins. Such an assay is also appropriate for “targeted discovery” experiments where the goal is to quantitate a large number (up to hundreds) of known proteins in a complex sample.

One approach for this strategy is the use of tandem mass spectrometry to monitor a unique peptide (or peptides) for each protein of interest by a selected reaction monitoring (SRM) assay, or by simultaneous analysis of many peptides by a multiple selected reaction monitoring (mSRM) assay. This approach can be extended further to provide absolute quantitation of targeted proteins by incorporation of appropriate stable isotope-labeled peptides as internal standards.

While mSRM assays are sensitive for targeted peptides, in a complex matrix, such as human serum, analyte selectivity can become a major issue. It is often difficult to differentiate between the targeted peptide signal and matrix background, particularly when quantifying multiple very low abundance proteins. The unique high resolution SRM (h-SRM) capability of the TSQ Quantum Ultra helps to restate this problem and increase assay specificity.

In this presentation, we demonstrate the TSQ Quantum Ultra mass spectrometer’s unparalleled capability for highly sensitive and accurate multiple protein quantitation from human serum by using high resolution multiple selected reaction monitoring (h-mSRM). An h-mSRM assay was developed for detecting 53 targeted proteins in human serum by using both unit mass resolution and high resolution (0.2 FWHM) for the Q1 quadrupole. The sensitivity, reproducibility, dynamic range and overall performance advantages of h-mSRM assays were evaluated. Additionally, a specific h-mSRM assay was developed for detecting a known biomarker (IL-6) from human serum.

Goal

Develop a fast, robust method for accurate, quantitative analysis of many targeted proteins in complex mixtures by using high resolution triple quadrupole mass spectrometry on a TSQ Quantum Ultra instrument.

Experimental Conditions

Sample Preparation

Whole human serum and interleukin-6 were used. The serum was diluted 40 times with 6M Guanidine. One milliliter of the diluted human serum sample and 10 µg of IL-6 were reduced and S-carboxymethylated, exchanged into 100 mM ammonium bicarbonate buffer and enzymatically digested. The digested mixtures were dried with a SpeedVac device and reconstituted with 200 µL water containing 0.1% TFA.

Peptide Selection and mSRM Transition Selection

Figure 1 shows two basic approaches for peptide selection and mSRM transition design. If the targeted protein is detected in a previous LC/MS/MS experiments, the peptides which a) had been detected repeatedly from these experiments, b) were unique for one single protein, c) contain no Cys, Met or other commonly modified residues and d) have proper mass range (600–2000 MW) were selected for the mSRM assays. Usually, multiple fragment ions for each selected peptide were used to maximize specificity.

If no HPLC/MS/MS data was available for the targeted protein, an SRM predictor tool (P3 Predictor) was used to predict candidate peptides and multiple fragment ions for SRM assay design (Figure 2). P3 Predictor takes amino acid sequences of targeted proteins of interest and performs in silico digestion. Peptides which contain no Cys, Met, His, N_xS(T) modification, R-P or K-P, and meet user-defined peptide length criteria will be listed as candidate peptides. A user simply selects one or multiple candidate peptides from the list and P3 Predictor will predict Q1 and Q3 SRM transitions automatically with proper collision energies and add them to an output csv file which the Quantum Ultra instrument can accept directly. In our experiments, for the 53 major serum proteins, 103 SRM transitions (Table 1) were used based on previous work by Anderson et al.⁽¹⁾ Six SRM transitions (560.82/616.38, 560.82/731.41, 560.82/844.49, 663.36/698.42, 663.36/812.46, 663.36/1012.54) were used for interleukin-6 based on the P3 Predictor tool. The collision energies were assigned by using the standard formula of $CE = 0.034 \times m/z + 3.314$.




Figure 1. Basic approaches for peptide selection and multiple SRM assay design

Paste Protein Sequence

MNSFSTSAFGPVAFSLGILLVLPAAFPAPVPPGEDSKDVAAPHRQPLTSSERID
 KQIRYILDGISALRKETCNKSNMCESSKEALAENNLNLPKMAEKDGCFCQSGFNE
 ETCLVKIITGLLEFVYLEYLQNRFSSEEQARAVQMSTKVLQFLQKKAKNLDAI
 TTPDPTTNASLLTKLQAQNLQLQDMTTHILRSFKELQSSLRALRQM

Peptide Features

Min Peptide Length: 7
 Max Peptide Length: 17

☒ Use Monoisotopic Precursor Ions
☒ Use Monoisotopic Product Ions

Exclude Peptides Containing:

☒ Cys
☒ Met
☒ His
☒ N-X/S/T motif
☒ R-P or K-P

Protein Features

☐ Exclude Potential Ragged Ends
 Eliminate first: 25 AAs
☐ Don't Digest

Check Library: None

Collision Energy Prediction
 Slope: 0.034
 Intercept: 3.314

Digest Protein
Reset

Peptide Sequences: 7

FFAPVPPGEDSK
 QPLTSSER
 YILDGISALR
 EALAENNLNLPK
 FESSEEQAR
 YLDGLQK
 EFLQSSLR

Precursor Info

Peptide: EFLQSSLR
 Monoisotopic Mass: M = 979.513
 M+H = 979.521
 [M+2H]²⁺ = 490.265
 [M+3H]³⁺ = 327.179
 Average Mass: M = 979.101
 M+H = 980.109
 [M+2H]²⁺ = 490.559
 [M+3H]³⁺ = 327.375
 Collision Energy [x1] = 36.6
 Collision Energy [x2] = 20
 Collision Energy [x3] = 14.4
 Hydrophobicity Retention Factor = 25.4

Product Ion Info

Seq #	B-Ion	#	Y-Ion
1	130.05	8	979.52
2	277.12	7	850.48
3	390.2	6	703.41
4	518.26	5	690.33
5	605.29	4	462.27
6	692.33	3	375.24
7	805.41	2	288.2
8	961.51	1	175.12

Predicted SRM Transitions

490.26	590.33	20
460.26	703.41	60
459.24	579.27	18.9
459.24	703.41	18.9
490.26	590.33	20
490.26	703.41	20
494.01	535.32	20.1
494.01	563.30	20.1
541.74	632.3	21.7
541.74	719.33	21.7
560.02	616.30	22.4
560.02	731.41	22.4
620.01	729.34	24.4
620.01	925.46	24.4
663.36	690.42	25.9
663.36	812.46	25.9

Add
Remove
Output CSV

Figure 2: P3 predictor—an multiple SRM transition prediction tool (by Michael J. MacCoss et. al from University of Washington)

LC Separation and MS Analysis

HPLC

A PicoFrit™ C18 column (75 µm×100 mm) was used for peptide separation. A Surveyor™ MS Pump was used to produce and deliver a solvent gradient (A:0.1%FA/2%ACN/98%H₂O, B:0.1%FA/100%ACN) to the column by means of a flow splitter. The post splitter flowrate was 300 nL/min. The linear ramp was from 2% B to 50% B in 85 min. Samples were loaded directly onto the column by a Micro AS autosampler after the flow splitter. The sample loading rate was 5 µL/min and loading time was 15 min.

MS

TSQ Quantum Ultra with Ion Max™ source equipped with a nanoflow column adapter (New Objectives) was used.

For SRM set up:

SRM set up 1: Q1, 0.7 FWHM; Q3, 0.7 FWHM
 SRM set up 2: Q1, 0.4 FWHM; Q3, 0.7 FWHM
 SRM set up 3: (h_SRM): Q1, 0.2 FWHM; Q3, 0.7 FWHM
 Q2: 1.5 mTorr (Ar)
 Scan width: 0.002 m/z
 Scan time: 20ms and 2 ms

For SRM-triggered MS/MS set up:

Scan Event 1
 Q1 and Q3: 0.7 FWHM; Q2: 1.5 mTorr;
 Scan width: 0.002 m/z; Scan time: 20 ms

Scan Event 2

DD precursor mass from Scan Event 1; Q1, 0.7 FWHM; signal threshold 30,000 counts
 Q2: 1.5 mTorr, CE: 0.034 × precursor mass m/z + 3.134;
 Dynamic Exclusion™ settings: repeat, 1; duration, 30s; exclusion time, 30s; exclusion list size, 50.

Protein	A	B	C	D	E	F	G	H	I	J	K	L	M
	Peptide Sequence	Q1	Q3	Run 1	Run 2	Run 3	Run 4	Run 5	Run 6	Ave	Std	CV (%)	
Area of total fragment ions													
Afamin	DADPDFFAK	563.8	826.4	460361	343876	453639	375549	358701	449877	407001	63231	13.1	
Alpha-1-acid glycoprotein 1	NWGLSVYADK-PETTK	563.8	940.4	3363991	3149662	4087824	3888070	3181300	3630443	3550348	386063	10.8	
Alpha-1-antichymotrypsin	EIGELYLPK	570.3	1052.5	3670363	3212687	3963132	3359817	3064500	3544359	3452476	316139	9.2	
Alpha-1B-glycoprotein	LETDFQLFK	570.3	633.4	3815348	3411667	4013810	3945637	3065305	3791761	3023935	219976	5.8	
Alpha-2-antiplasmin	LGNDEPGGGTALK	566.8	771.4	1063920	804714	1221469	1097863	812108	1013693	1002295	165111	16.5	
Alpha-1-antitrypsin	DTEEDFHYDGVTVTK	570.3	790.4	2103352	2725101	1671968	1594758	1806790	182299	1680711	841187	50	
Alpha-2-macroglobulin	LLIYAVLTGDVIGDSAK	570.3	1172.6	9115043	7565135	6513412	9040774	8044853	7233501	7918786	1027185	13	
Angiotensinogen	PKDPTFIPAPIQAK	508.3	556.4	157986	222990	226592	226163	163494	209165	201065	31923	15.9	
Antithrombin-III	DDLVSDFAHK	437.2	803.4	208027	274147	172331	143431	167764	147782	185580	49084	25.4	
Apolipoprotein A-I	ATEHLSTLSEK	405.9	664.4	8310093	7404901	8768816	7389413	8732920	8181190	8129556	611263	7.5	
Apolipoprotein A-II precursor	SPELQAEAK	486.8	546.4	14652243	13161798	16449486	15604447	15516303	15091455	15084289	1111632	7.4	
Apolipoprotein A-IV	SLAPYAGDTQEK	675.8	982.4	882356	1040338	943183	1086796	1066425	978521	999603	78869	7.9	
Apolipoprotein B-100	FPEYDVLTK	524.3	803.5	815385	751591	822726	844726	590265	666462	748524	101016	13.5	
Apolipoprotein C-I liprotein	TPDVSSALDK	516.8	719.3	3775528	3585466	4138996	3724512	3760057	4077021	3841930	217222	5.7	
Apolipoprotein C-II liprotein	STAAASTYTGIFDQVLSVLK	760.4	1149.7	6767	13319	47439	20992	46405	37616	26760	17437	60.6	
Apolipoprotein C-III liprotein	DALSSVQESQVAGQAR	750.4	1002.6	1508129	1271980	1643745	1482611	1323109	1294646	1420737	147500	10.4	
Apolipoprotein E	LGPLVEGGR	484.8	588.3	1219946	1091977	1272193	1028347	1144104	1103460	1143338	89358	7.8	
Beta-2-glycoprotein I	ATVYVGGR	511.8	751.4	9689076	8814814	10796000	9447598	10202994	10095120	9874267	679873	6.9	
C4b-binding protein alpha chain	LSLIEIQLELQR	511.8	652.3	1350409	1340156	963933	1775467	1643076	1614620	1447944	292765	20.2	
Ceruloplasmin	EYTDASFNR	735.9	815.5	1856971	1680288	1787627	2014591	2157404	1910802	1901281	168870	8.9	
Clusterin	LFDSDPITVTPVEYSR	602.3	824.3	2395112	2520417	2362234	2690174	2635699	2777964	2566915	161392	6.3	
Coagulation factor V	DPSPDLLLK	937.5	686.4	2957516	3484642	4611308	2649725	3077363	3572073	3392105	687963	20.3	
Coagulation factor XIIIa heavy chain	VVGLVALR	665.8	898.6	588570	518838	828562	898422	898715	892279	733964	106135	14.5	
Complement C3	TGLQEVVK	442.3	686.4	4566881	4737491	5720918	5120795	5411222	5561088	5186389	481795	8.9	
Complement C4 gamma chain	ITQVLHFK	501.8	603.3	111095	112564	100689	159819	133237	128846	124325	21141	17	
Complement C4 beta chain	VGDTLNLR	362.9	843.5	2302932	2038312	1821734	2460076	1632647	2100302	2059334	303503	14.7	
Complement factor C9	AIEDYNEFSVR	557.8	629.4	1271.6	185702	212040	174997	161374	230387	205447	194991	25570	13.1
Complement factor B	EELLPAQDIK	728.5	1027.5	2011175	2312117	2128099	2011175	2447949	2128099	2173102	173987	8	
Complement factor H	SPDVNSPISQK	578.4	784.5	330940	425920	267351	267545	427823	267351	331155	78120	23.6	
Fibrinogen alpha chain	GSESGIFNTK	671.4	830.4	110715	94678	137387	116792	138690	138164	122738	18299	14.9	
Fibrinogen beta chain	QGFGNVATNTDGK	570.8	887.5	64139	61523	76351	87158	91996	88959	78354	13156	16.8	
Fibrinogen gamma chain	DTVQIHDTGK	654.8	705.4	980969	1127408	1211679	1290416	1105767	1122725	1139811	104564	9.2	
Fibronectin	DLOFVEVDVK	409.5	533.3	148897	79771	95323	118285	113737	124585	113433	23945	21.1	
Gelsolin, isoform 1	TGAQELLR	647.3	680.4	422213	430807	504372	465432	559462	569154	491907	63210	12.9	
Lactoglobulin beta chain	VGYSVGWR	444.3	786.5	36340479	30692289	43475501	35250502	44157166	46000993	40652022	4402765	11	
Leiomopexin	NFPSPVDAAFR	490.6	661.3	19718474	17178303	21208304	18586954	17250682	21933002	19312620	1997910	10.3	
Leiparin cofactor II	TLEAQLTPR	610.8	959.6	1176390	1501771	1764117	1165217	1594067	1725999	1487927	262912	17.7	
Lysidine-rich glycoprotein	DSPVLIDFFEDTER	514.8	814.4	1249897	1304182	1606963	1714889	1660153	1811153	1557873	228422	14.7	
Inter-alpha-trypsin inhibitor heavy chain	AAISGENAGLVR	841.9	1171.5	2118520	2104894	2392357	2244207	2302923	2245588	2234748	109605	4.9	
Inter-alpha-trypsin inhibitor light	AFIQLWAFDAVK	941.9	1058.4	1316403	1797675	1834854	1347052	1395863	1588799	1546743	229564	14.8	
Kininogen	TVGSDTFYSFK	704.9	940.5	2601441	2998187	3752582	2872619	2383236	3727717	3055964	571426	18.7	
Le-selectin	AEIEYLEK	626.3	994.5	77482	105168	91786	96749	102618	97221	95171	9864	10.4	
Plasma retinol-binding protein precursor	YWGVSFLQK	497.8	881.3	2178485	1964456	2480949	2284090	2047641	2444580	2233367	208950	9.4	
Plasminogen	LFLEPTR	559.8	849.5	1964456	2480949	2284090	2047641	2444580	2233367	208950	9.4		
Prothrombin	ETAASLLQAGYK	438.3	615.4	2221967	1958030	2432750	2097034	1719680	2008048	2072918	242677	11.7	
Serum albumin	LVNEVFPAK	438.3	502.3	1543698	1348273	1748347	1559187	1802953	1570716	1562196	128452	8.2	
Serum amyloid P-component	VGEYSLVGR	626.3	879.5	580114007	506210531	661792006	512338571	503622280	583128349	5.58E+08	62637363	11.2	
Transferrin	EDPQTIFYAVAVVK	575.4	937.4	504502	508405	563431	552061	554938	548342	37744	6.9		
Transferrin	EDPQTIFYAVAVVK	575.4	694.4	504502	508405	563431	552061	554938	548342	37744	6.9		
Transferrin	EDPQTIFYAVAVVK	575.4	937.4	504502	508405	563431	552061	554938	548342	37744	6.9		
Vitamin D-binding protein	THLPEVFLSK	575.4	871.5	1797675	1834854	1347052	1395863	1588799	1546743	229564	14.8		
Vitamin K-dependent protein C	WELDLDIK	815.4	932.5	168787	145466	258455	158618	200416	189570	189570	41128	21.7	
Vitronectin	FEDGVLDPDYPR	516.3	603.3	1924890	1666834	1919119	1972180	1809208	2005032	1882877	124962	6.6	
Zinc-alpha-2-glycoprotein	EIPAWVFPDPAAGITK	711.9	1031.5	2062315	2009293	2619535	2236890	1759158	2666391	2225430	358207	16.1	

Table 1. Multiple SRM assays and reproducibility for 53 major serum proteins with TSQ Quantum Ultra

Results

Quantitative results for 53 targeted major serum proteins

A total of 103 unique SRM transitions were chosen as proteotypic peptides for the 53 targeted proteins of interest. Using identical human serum samples, five nanoflow HPLC-MS/MS experiments were performed for each sample, with the same 103 SRMs monitored with different resolutions (0.2, 0.4 and 0.7 FWHM) and Q2 dwell times (2 and 20 ms). The high resolution SRM assays (Q1: 0.2 FWHM) gave the best results and clearly resolved targeted analyte transitions from interference peaks that were seen at lower Q1 resolutions. Figure 3 shows one example where high resolution was used to unambiguously detect peptide QGFGNVATNTDGK

(representing fibrinogen beta chain). Importantly, it should be noted that 0.2 FWHM data revealed the presence of significant matrix interference in 25% of the SRM transitions which were monitored. The 103 SRMs were detected with enough scans for reliable quantitation using both 2 ms and 20 ms scan times, although the 2 ms scan time gave twice as many scans with only a minor decrease in signal intensity (Figure 4). For narrow peak widths, such as those associated with uHPLC, 2 ms scan times will be required.

For testing this method's reproducibility, the same h-mSRM experiment at Q1: 0.7 FWHM and 20 ms scan time was repeated six times and the results were summarized in Table 1. Among the 53 targeted proteins,




Figure 3. SRM traces for peptide QGFGNVATNTDGK (fibrinogen beta chain) at different Q1 resolution settings, showing the power of high resolution SRM (h-SRM) for unambiguous detection of targeted peptides from human serum. Note the marked specificity of the Q1 0.2 FWHM trace relative to those at 0.4 and 0.7. This performance is achieved with very little reduction in signal intensity.




Figure 4: SRM traces for peptide LGPLVEQGR with different SRM dwell times

51 proteins produced acceptable quantitative data, while only two were not reliably observed. For the whole serum digests, CVs (n=6) were from 5–26% (50% of SRMs had CVs 10%). Proteins present at concentrations down to $\mu\text{g/mL}$ levels, such as L-selectin and fibronectin,⁽¹⁾ were reliably detected, yielding a dynamic range of greater than four orders of magnitude (from lowest peak areas of

$6\text{E}+04$ from fibrinogen beta chain to the highest peak areas of $6\text{E}+08$ from albumin peptide) in a single experiment. MS/MS spectra were acquired once the SRM intensity exceeded 30,000 counts and typically of good quality and showed rich y series and some b series fragment ions, permitting database searching with SEQUEST® in BioWorks™ 3.3 (Figure 5).




Figure 5: SRM transition for peptide ATEHSTLSEK (Apolipoprotein A-I) triggered the MS/MS spectrum shown below. Database searching with SEQUEST in BioWorks confirms that the peptide structure agrees with the predicted transition.




Figure 6: Detection limit and linear dynamic range of IL-6 spiked into human serum

Detection limit and linear dynamic range of IL-6, using h-mSRM assays developed with P3 Predictor

An IL-6 digest mixture was spiked into human serum to create a dilution series, and was run sequentially to evaluate the detection limit and linearity of the h-mSRM assay. Excellent analytical sensitivity and linearity were seen (Figure 6). The linear dynamic range for spiked IL-6 was over four orders of magnitude (0.02–400 fmol on column) and the limit of detection was 20 amol on column with a S/N of 22.

Conclusions

An h-mSRM- based assay for targeted proteins in human serum was developed on a TSQ Quantum Ultra triple quadrupole mass spectrometer. A total of 103 SRM transitions were monitored for the quantitation of 53 serum proteins. One exogenous protein, IL-6, was spiked to evaluate assay performance. Several unique features of the Quantum Ultra instrument contributed to the specificity, sensitivity, and robustness of this assay.

1) Q1 resolution of 0.2 FWHM dramatically reduced non-specific matrix interference from the serum background, dramatically improving assay specificity. At Q1 resolution of 0.4 or 0.7 FWHM, significant interference was seen in 25% of targeted transitions.

2) Analytical assay performance was excellent.

- %CV varied from 5-26%, with 50% of protein CVs < 10%.
- Peptide response was linear over four orders of magnitude.
- Sensitivity was excellent, with the ability to detect proteins present at µg/mL levels and to detect IL-6 at levels as low as 20 amol on column.
- SRM-triggered MS/MS spectra were of good quality and in most cases sufficient to permit confirmation of peptide ID by database searching with SEQUEST.

References

- ⁽¹⁾ Leigh Anderson and Christie L. Hunter (2006) Quantitative Mass Spectrometric Multiple Reaction Monitoring Assays for Major Plasma Proteins. *Mol. Cell. Proteomics* 5.4, 573-588.
- ⁽²⁾ Michael J. MacCoss et. al. (2006) Private communication.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

Australia
+61 2 8844 9500
Austria
+43 1 333 50340
Belgium
+32 2 482 30 30
Canada
+1 800 532 4752
China
+86 10 5850 3588
Denmark
+45 70 23 62 60
France
+33 1 60 92 48 00
Germany
+49 6103 408 1014
India
+91 22 6742 9434
Italy
+39 02 950 591
Japan
+81 45 453 9100
Latin America
+1 608 276 5659
Netherlands
+31 76 587 98 88
South Africa
+27 11 570 1840
Spain
+34 91 657 4930
Sweden/Norway/Finland
+46 8 556 468 00
Switzerland
+41 61 48784 00
UK
+44 1442 233555
USA
+1 800 532 4752

www.thermo.com

©2007 Thermo Fisher Scientific Inc. All rights reserved. SEQUEST is a registered trademark of the University of Washington. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



Thermo Finnigan LLC,
San Jose, CA USA is ISO Certified.
AN62260_E 01/07S

View additional Thermo Scientific LC/MS application notes at: www.thermo.com/appnotes