

Highlighting innovative design features
and useful application information for
Thermo Scientific Centrifuges and Rotors

Thermo
S C I E N T I F I C

smart notes

► design and innovation



LARGE CAPACITY FLOOR MODEL CENTRIFUGES

SMART NOTE 1

Q A

How can blood banks eliminate potential variables to ensure consistent product yield and reproducibility?

Traditional centrifuge run settings – speed and time – do not take into consideration instrument and environmental factors, such as centrifuge load, instrument age, or altitude, and as a result, reproducibility can be compromised. The Thermo Scientific Accumulated Centrifugal Effect (ACE) integrator function automatically compensates for variations by ensuring that separating g-forces are calculated and adjusted during run time for protocol consistency – across runs, labs and even sites.

Run-to-run repeatability and product consistency are important objectives of every blood component facility. In fact, the first centrifugation step for platelet production is critical to the final result. For example, if the run is too long, platelets separate into red blood cells and buffy coat. Conversely, if the run is too short, red and white cells are not separated completely from platelet-rich plasma. The ACE™ integrator function on Thermo Scientific large capacity centrifuges eliminates potential variables to ensure consistent run performance.



Run-to-run repeatability and product consistency are critical for every blood component facility.

How ACE Integrator Function Works

Centrifugation runs involve speed and time, along with variables influencing instrument performance, such as full or partial rotor load, voltage fluctuations, loss of instrument calibration, and environmental factors including altitude and extreme ambient temperatures. The ACE integrator function calculates the g-force experienced during the run in increments of speed over time to give a value representing the overall separating g-force. This value can be substituted for the "TIME" setting, therefore duplicating the overall separating g-force for every run.

Best Practices for Product Consistency

Obtaining a consistent product requires understanding and controlling process variables. By eliminating variables inherent in centrifugation, such as rotor load, the ACE integrator function allows best practices to be established and maintained from operator-to-operator, run-to-run, instrument-to-instrument, and site-to-site. The ACE integrator function is a useful and powerful feature of Thermo Scientific large capacity centrifuges to ensure process control and ultimately sample protection.

► Summary

The ACE integrator function on Thermo Scientific large capacity centrifuges eliminates potential variables to ensure consistent product yield and reproducibility.

Results Without the ACE Integrator Function



Figure 1: In a typical first centrifugation step, a two-bag rotor load attains set speed faster than a six-bag load. Since both loads will time out at the set time of 3:30 minutes, different g-forces are achieved during the run. By using the ACE integrator function, the time for 2 bags would be changed to 3:00 minutes to obtain the same overall g-force for both loads.

Results With the ACE Integrator Function




Figure 2: With an ACE value and speed set at the start of a run, times were adjusted to achieve the same overall g-force regardless of the rotor load.

Experience our large capacity centrifuges now at www.thermoscientific.com/centrifugevideos

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America: USA/Canada +1 866 984 3766 (866-9-THERMO)

www.thermoscientific.com/centrifuge

Europe: Austria +43 1 801 40 0, Belgium +32 53 73 42 41, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 95059 448, Netherlands +31 76 579 55 55, Nordic/Baltic/CIS countries +358 9 329 10200, Russia +7 812 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: Australia +61 39757 4300, China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220, New Zealand +64 9 980 6700, Other Asian countries +852 2885 4613 **Countries not listed:** +49 6184 90 6940