

Isolation of RNA Using the Thermo Scientific Sorvall MTX Micro-Ultracentrifuge and Rotors

KEY WORDS

- RNA Isolation
- Cell Culture
- Sorvall MTX Micro-Ultracentrifuge
- Sorvall MX Micro-Ultracentrifuge
- S55-S Swinging Bucket Rotor
- S150-AT Fixed Angle Rotor

Introduction

Many protocols exist for the extraction and purification of total RNA from cells or tissue. The key to obtaining high quality RNA is the rapid lysis of cells, the thorough denaturation of endogenous nucleases, and ultimately, the efficient removal of these nucleases.

Early methods of RNA purification relied on the use of guanidinium thiocyanate (GTC) for the rapid denaturation of cellular proteins and extraction of RNA by ultracentrifugation through a cesium chloride (CsCl) cushion or a single step procedure that uses a mixture of GTC and phenol-chloroform. Recent techniques and commercially available kits have been developed around these methods that allow for a more rapid process for RNA isolation.

Using microcentrifuges, these methods have been found to decrease the time required for RNA isolation from a few hours to less than an hour. This brief will provide protocols for the traditional methods of RNA isolation using fixed-angle and swinging bucket rotors in a Thermo Scientific Sorvall MTX Micro-Ultracentrifuge, and will also identify commercially available kits that are compatible with Thermo Scientific microcentrifuges.

Procedures

Protocol 1:

Rapid Isolation of RNA from Cultured Cells Using the Thermo Scientific S55-S Swinging Bucket Rotor and Sorvall® MTX Micro-Ultracentrifuge^{3,4}

This procedure describes the use of GTC and ultracentrifugation through a CsCl cushion with a S55-S swinging bucket rotor to isolate RNA from a cell culture. This procedure has been developed for working with small samples and offers a high yield (>90%) of RNA in a short period of time (approximately 3 hours).

Lysis

1. Remove media from cell culture plate.
2. Add 600-1000 μL of GTC solution (4.0 M guanidine thiocyanate, 0.03 M sodium acetate, 1% β -mercaptoethanol) to each well of the culture plate to lyse cells.
3. Lift cells with a cell scraper and homogenize the slurry by aspirating 4-5 times through a 22 gauge needle.



Thermo Scientific S150-AT Rotor

Isolation

1. Pre-load 2.2 mL polyallomer (PA) tubes (PN 45240) with 550 μL of 5.7 M CsCl and mark solution level with a waterproof pen (Figure 1).
2. Layer samples onto the CsCl cushion and fill tubes to within 3 mm of the top of the tube with the GTC solution.
3. Place tubes in a S55-S swinging bucket rotor in a Sorvall MTX Micro-Ultracentrifuge with the following parameters: 55,000 RPM (258,826 $\times g$) for 3.0 – 3.2 hours at 22 °C, a Thermo Scientific Sorvall MX Micro-Ultracentrifuge may alternatively be used. *Note: Non-precipitating solutions up to 1.7 g/mL can be run in the S55-S rotor without a reduction in rotor speed. Consult the rotor manual for instructions on rotor deration and how to reduce the chance of CsCl precipitation.*
4. Following ultracentrifugation, use a sterile Pasteur pipette to remove the top layer of solution down to the interface reference mark (Figure 2).
5. Using a second sterile pipette, remove the remaining solution down to approximately 250 μL .
6. Carefully cut off the top of each tube just below the reference mark, pour off the remaining solution, and blot the tubes dry on a clean paper towel.



Thermo Scientific Sorvall MTX Micro-Ultracentrifuge

Recovery

1. Rinse the RNA pellet X 2 with 75 μ L of 70% isopropanol and blot the tubes dry on a clean paper towel.
2. Add 100 μ L of DEPC-treated water to each tube and freeze (-70 °C) for a minimum of 45-60 minutes. *Note: DEPC is a suspected carcinogen; follow manufacturers recommended handling procedures.*
3. Resuspend RNA by pipetting up and down several times.
4. Transfer the solution to a sterile RNase-free 1.5 mL microcentrifuge tube.
5. Wash the cut ultracentrifuge tube with a second 100 μ L aliquot of 0.25 volumes of 7.5 M ammonium acetate and 2.5 volumes of 100% ethanol.
6. Collect precipitated RNA by centrifugation at >10,000 x g for 20 minutes at 4 °C in a Thermo Scientific refrigerated microcentrifuge.
7. Wash the pellet with cold 70% ethanol and vacuum dry.
8. Resuspend the RNA pellet in a small volume of DEPC-treated water.

Protocol 2:

Separation of Total RNA from a Homogenated Cell Line with the S150-AT Fixed Angle Rotor in Less than 2 Hours^{4,5}

This procedure describes the separation of total RNA from a homogenate cell line using the Sorvall MTX Micro-Ultracentrifuge and the S150-AT fixed angle rotor, which has a maximum speed of 150,000 rpm. This procedure can separate total RNA in only 1.75 hours.

Isolation

1. Pellet cells (5-10 \times 10⁶ cells/tube).
2. Add 1 – 1.5 mL of 4 M GTC and vortex to the lyse cells.
3. Layer the samples onto a 0.8 mL 5.3 M CsCl solution that has been pre-loaded in a 2.0 mL polyallomer (PA) Re-Seal™ tube (PN 45246).
4. Place tubes in the S150-AT fixed angle rotor and centrifuge in a Sorvall MTX Micro-Ultracentrifuge with the following parameters:




Figure 1




Figure 2

Figure 1. Polyallomer tube containing the sample (GTC solution) on top of the CsCl cushion.

Figure 2. Withdrawal of the supernatant following centrifugation. The figure illustrates the depth of the pipette tips and the cut position relative to the reference mark placed on the tube delineating the interface between the sample and the CsCl cushion.

5. Dissolve the RNA pellet in DEPC-treated water and perform ethanol precipitation. *Note: Separation of total RNA may not be completed due to crystallization during centrifugation when 5.7M CsCl is used.*

Conclusion

The isolation and purification of RNA can be accomplished using multiple methods. This brief describes traditional approaches for RNA purification using ultracentrifugation with swinging bucket and fixed-angle rotors. A swinging bucket rotor is often used because the RNA is collected at the bottom of the tube. Although centrifugation in a vertical or fixed angle rotor requires less time for RNA purification, it may result in materials collecting or coming in contact with the sides of the tube, thus potentially contaminating the RNA. Alternatively, commercially available kits (Promega®, Ambion®, and Qiagen™) have been developed that are compatible with Thermo Scientific microcentrifuges and may offer enhanced time savings, while increasing RNA yield and purity.

References

1. Chirgwin, J.J., A.E. Przybyla, R.J. MacDonald, and W.J. Rutter (1979). Isolation of Biologically Active Ribonucleic Acid From Sources Enriched in Ribonuclease. *Biochemistry* 18:5294.
2. Glisin, V., R. Crkvenjakov, and C. Byrus (1973). Ribonucleic Acid Isolated by Cesium Chloride Centrifugation. *Biochemistry* 13:2633.
3. Norton, J. Rapid Isolation of RNA from Cultured Cells Using the S55-S Swinging Bucket Rotor and Discovery Microultracentrifuge. Application Brief, S00266.
4. Sambrook, J. and D. Russell (2001). *Molecular Cloning: a laboratory manual*, 3rd ed. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, New York.
5. Separation of total RNAs from a homogenated cell line with the S150-AT fixed angle rotor in less than 2 hours. Diverse Applications for Discovery M120 & M150 Micro-ultracentrifuges. Application Brief, S00242.

In addition to these offices, Thermo Fisher Scientific maintains a network of representative organizations throughout the world.

North America:

USA / Canada
+1 866 984 3766
(866-9-THERMO)

Europe:

Austria
+43 1 801 40 0

Belgium
+32 53 73 42 41

France
+33 2 2803 2180

Germany national toll free
08001-536 376

Germany international
+49 6184 90 6940

Italy
+39 02 95059 448

Netherlands
+31 76 579 55 55

Nordic/Baltic/CIS countries
+358 9 329 10200

Russia
+7 812 703 42 15

Spain / Portugal
+34 93 223 09 18

Switzerland
+41 44 454 12 12

UK / Ireland
+44 870 609 9203

Australia
+61 39757 4300

Asia:
China
+86 21 6865 4588 or
+86 10 8419 3588

India toll free
1800 22 8374

India
+91 22 6716 2200

Japan
+81 45 453 9220

New Zealand
+64 9 980 6700

Other Asian countries
+852 2885 4613

Countries not listed:
+49 6184 90 6940

www.thermoscientific.com/centrifuge

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Re-Seal is a trademark of Seton Scientific. Promega is a registered trademark of Promega Corporation. Ambion is a registered trademark of Ambion, Inc. Qiagen is a trademark of Qiagen. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.