

Utilization and Optimization of Cryopreserved Human Hepatocytes as a Model to Assess CYP450 Inhibition

LaHoma Easterwood, Helen Rollison¹, Lee Robinson, Ben Fontenot, Jenny Smith, Mike Shaw, and Jason Wright

Cell Systems Division • Life Technologies • 1624 Headway Circle • Austin, Texas 78754 • USA and (1) Clinical Pharmacology and DMPK, Alderley Park, AstraZeneca, Cheshire SK10 4TG, UK




Introduction

Pharmacokinetic drug-drug interactions (DDIs) can occur when one drug alters the metabolism of a co-administered drug. To gain a better understanding of a drug's potential to cause pharmacokinetic-based DDIs via inhibition of drug metabolizing enzymes, a battery of *in vitro* screens can be employed to determine the extent of CYP inhibition. These assays typically use one or more of the following tools: recombinantly expressed CYPs (rCYPs) or human liver microsomes (HLM). Human primary hepatocytes, which provide an *in vitro* environment which more closely resembles that of the human liver by providing the full complement of Phase I and Phase II xenobiotic-metabolizing enzymes, have recently been suggested as a model to determine both reversible and time-dependent inhibition of CYPs. For this study we utilized cryopreserved human hepatocytes in suspension as a tool to assess reversible and time-dependent inhibition of CYP enzymes by isoform-specific inhibitors. The inhibition potential of these compounds for the major drug metabolizing CYP enzymes (CYP1A2, phenacetin *O*-dealkylation; CYP2B6, bupropion hydroxylation; CYP2C8, paclitaxel 6 α -hydroxylation; CYP2C9, tolbutamide hydroxylation; CYP2C19, (S)-mephenytoin 4'-hydroxylation; CYP2D6, dextromethorphan demethylation; 6-hydroxylation; CYP3A4, testosterone 6 β -hydroxylation and midazolam 1'-hydroxylation) was assessed by measuring IC_{50} values for direct inhibition as well as K_i and k_{inact} values for time-dependent inhibition using known isoform-specific positive control inhibitors. Prior to determining the inhibition potential, the reactions conditions were optimized by assessing the linearity of time and cell density dependence and the K_m and V_{max} parameters for each of the isoforms were determined. These resulting kinetic parameters obtained are comparable to those determined for human liver microsomes and demonstrate that cryopreserved hepatocytes can be used for assessing both direct and time-dependent CYP inhibition potential by drug candidates.



Figure 1 – Linearity assessment



Linearity of 1-hydroxymidazolam formation as a representative dataset for linearity assessment. Incubations for CYP3A4 (midazolam) were performed for 120 minutes at 37 °C


Linearity of metabolite formation from each of the probe substrates was assessed over a range of 10 to 120 minutes incubation and at cell concentrations of 0.1, 0.25 and 0.5 million cells per mL. Each probe substrate was incubated at 2 to 5 times the known K_m concentration in human liver microsomes. Incubations were started by the addition of hepatocytes and terminated by the addition of the appropriate stop solution containing internal standard. Metabolite formation was analysed using the HPLC-MS/MS methods from the CellDirect validated human liver microsomal assays. Linear conditions were determined for each CYP isoform and all the HPLC-MS/MS methods were judged suitable for use.

Table 1 – Conditions used for incubations in cryopreserved human hepatocytes

| Human CYP | Substrate | Substrate concentration (μ M) | Cell concentration (million cells per mL) | Incubation time (minutes) |
|-----------|------------------|------------------------------------|---|---------------------------|
| CYP1A2 | Phenacetin | 50 | 0.5 | 120 |
| CYP2B6 | Bupropion | 60 | 0.1 | 120 |
| CYP2C8 | Paclitaxel | 5 | 0.25 | 120 |
| CYP2C9 | Diclofenac | 5 | 0.25 | 30 |
| CYP2C19 | (S)-mephenytoin | 50 | 0.25 | 120 |
| CYP2D6 | Dextromethorphan | 5 | 0.1 | 120 |
| CYP3A4/5 | Midazolam | 5 | 0.1 | 120 |
| CYP3A4/5 | Testosterone | 50 | 0.1 | 120 |

Incubations were performed under linear conditions at the approximate K_m concentration for each probe substrate. CellDirect/Life Technologies cryopreserved human hepatocytes designated as Hu800, Hu8058, and Hu8066 were used in all incubations.

Figure 2 – Kinetics assessment



Michaelis-Menten plot of 1-hydroxymidazolam formation as a representative dataset for the kinetics assessment. Incubations for CYP3A4 (midazolam) were performed for 120 minutes at 37 °C

Incubations for all isoforms assessed were performed under linear conditions for each probe substrate with hepatocyte concentrations of 0.1 – 0.5 million cells per mL and incubation times of 30 or 120 minutes. Variable concentrations of each probe substrate were used (Five (5) or six (6), depending on assay) and samples were terminated and analysed as described for the linearity assessment. Systat was used to perform the non-linear regression of the rates of metabolite formation at each substrate concentration. The kinetic parameters, K_m and V_{max} , were generated from the fit of the dataset.


Table 2 – Comparison of kinetic parameters in cryopreserved human hepatocytes versus human liver microsomes

| Human CYP | Substrate | K_m value in hepatocytes (μ M) | K_m value in HLM (μ M) |
|-----------|------------------|---------------------------------------|-------------------------------|
| CYP1A2 | Phenacetin | 38 | 40 |
| CYP2B6 | Bupropion | 11 | 56 |
| CYP2C8 | Paclitaxel | 5.1 | 9.3 |
| CYP2C9 | Diclofenac | 40 | 5.3 |
| CYP2C19 | (S)-mephenytoin | 13 | 33 |
| CYP2D6 | Dextromethorphan | 2.3 | 3.7 |
| CYP3A4/5 | Midazolam | 2.6 | 2.8 |
| CYP3A4/5 | Testosterone | 78 | 60 |

Table 3 – Percentage inhibition in cryopreserved hepatocytes when co-incubated with positive control inhibitors

| Human CYP | Chemical inhibitor | Inhibitor concentration (μ M) | % Inhibition |
|-----------|--------------------|------------------------------------|--------------|
| CYP1A2 | Furafylline | 5 | 42.2 |
| CYP2B6 | ThioTEPA | 20 | 47.0 |
| CYP2C8 | Quercetin | 10 | 11.9 |
| CYP2C9 | Sulfaphenazole | 6.3 | 62.7 |
| CYP2C19 | Ticlopidine | 1 | 73.2 |
| CYP2D6 | Quinidine | 0.4 | 41.4 |
| CYP3A4/5 | Ketoconazole | 0.1 | 58.4 |
| CYP3A4/5 | Ketoconazole | 0.1 | 39.8 |

Figure 3 – Inhibition of CYP3A4/5 in cryopreserved hepatocytes by mifepristone



Inhibition parameters of $K_i = 23.1 \mu$ M and $k_{inact} = 0.101 \text{ min}^{-1}$ were determined for the inhibition of CYP3A4/5 by mifepristone.

Table 4 – Percentage inhibition in cryopreserved hepatocytes when pre-incubated with positive control inhibitors

| Human CYP | Chemical inhibitor | Inhibitor conc (μ M) | % Inhibition |
|-----------|--------------------|---------------------------|--------------|
| CYP1A2 | Furafylline | 1 | 70.7 |
| CYP2B6 | ThioTEPA | 30 | 64.5 |
| CYP2C8 | Phenelzine | 100 | 30.9 |
| CYP2C9 | Tienilic acid | 3 | 41.7 |
| CYP2C19 | Ticlopidine | 0.5 | 60.6 |
| CYP2D6 | MDMA | 10 | 44.6 |
| CYP3A4/5 | Mifepristone | 10 | 64.9 |
| CYP3A4/5 | Mifepristone | 10 | 54.7 |

Results and Conclusions

- Hepatocytes can be a useful tool to determine both reversible and time-dependent inhibition of CYPs.
- The use of hepatocytes provides an *in vitro* environment which more closely resembles that of the human liver by providing the full complement of Phase I and Phase II xenobiotic-metabolizing enzymes.
- Assays for the identification of reversible and time-dependent inhibition of CYP 1A2, 2B6, 2C8, 2C9, 2C19, 2D6 and 3A4/5 were successfully developed and validated.
- The kinetic parameters determined for CYP 1A2, 2B6, 2C8, 2C9, 2C19, 2D6 and 3A4/5 in cryopreserved human hepatocytes were comparable to those determined for human liver microsomes.
- The utility of using hepatocytes was shown through the development and validation of a time-dependent assay for CYP3A4/5 using mifepristone as the positive control.