

Liquid Chromatography Optimization Glycopeptides Analysis by Electron Transfer Dissociation

Terry Zhang, Rosa Viner and Vlad Zabrouskov ; Thermo Fisher Scientific, San Jose, CA

Overview

Purpose: Optimize the liquid chromatography for glycopeptides analysis by Electron Transfer Dissociation (ETD).

Methods: Glycoprotein digests were analyzed by nanoLC ESI using different types of HPLC columns and a Thermo Scientific LTQ XL™ mass spectrometer equipped with an ETD source option.

Results: Graphitic carbon column demonstrated excellent capabilities for glycopeptides analysis by mass spectrometry without any enrichment.

Introduction

Glycosylation plays key roles in controlling various biological processes. It is one of the most widespread and complex post-translational modifications (PTMs) found on proteins and its characterization remains a great analytical challenge. LC-MS/MS is the most powerful and versatile technique for glycopeptide structure elucidation. However, commonly used collisional-induced dissociation (CID) has limitations on determining the modification site due to the labile nature of the glycan modification. Also sialic acid residue in glycans of glycoproteins add negative charges, and thus reduces the pI of the corresponding molecule. This additionally complicates the glycoprotein analysis by LC/MS. As a new dissociation technique, Electron Transfer Dissociation (ETD) preserves labile PTMs and provides a new and powerful tool that makes the identification of modification site possible. Since glycosylated proteins and the resulting peptides are often highly heterogeneous, the high quality liquid chromatography is critical for glycopeptides analysis.

In this study, two reasonably well characterized glycoproteins, bovine α₁-acid glycoprotein and bovine fetuin [1,2] were analyzed using nano LC MS/MS with ETD. Liquid chromatography separation conditions were systematically optimized for glycopeptides analysis prior mass spectrometry using different stationary-phase columns. Two reversed phase columns, C₈ and C₁₈, and a graphitic carbon column were evaluated.

Methods

Samples:

Reduced and alkylated enzymatic digests of bovine α₁-acid glycoprotein, bovine fetuin. Glycoproteins were purchased from Sigma.

LC/MS:	
HPLC System:	Surveyor™ MS Pump with a flow splitter
Columns:	Agilent® ZORBAX 300SB C8 column (75 μm x 5 cm) Microtech, C18 column (150 μm x 10cm) Thermo Scientific, Hypercarb™ column (75 μm x 5cm)
Mobil Phase:	A: Water, 0.1% formic acid; B: Acetonitrile, 0.1% formic acid
Gradient:	5-50% B in 30 minutes
Mass Spectrometer:	Thermo Scientific LTQ XL™ linear ion trap mass spectrometer with ETD and nano-ESI source
Spray Voltage:	2 kV
Capillary Temp:	160 °C
Capillary Voltage:	35 V
Tube Lens:	125 V
MSn Target:	1e4
Mass range :	50-2000 m/z or 100-4000 m/z
Anion Reagent:	Fluoranthene
Anion Reagent Isolation:	On
Anion Target:	2e5
Max Anion Injection Time:	50 ms
ETD Reaction Time:	100 ms
Mass Spectrometer	Thermo Scientific LTQ Orbitrap XL™
Mass range	400-2000 m/z, resolution 100000

Data Processing:

Data were processed using BioWorks™ 3.3.1 with SEQUEST®. Xtract program (Thermo Scientific) was used for deconvolution of multiply charged precursors.

Results

Figure 1 shows LC/MS analysis of bovine α₁-acid glycoprotein on C₁₈ and Hypercarb columns. The top panels are the base peak chromatograms and the bottom panels are the bi-antennary glycopeptide 91-99 MS profiles. One pmol of protein digest was injected into C₁₈ column versus 500 fmol on Hypercarb column resulting in a three times higher C₁₈ base peak intensity than the one from the Hypercarb column. However as shown in bottom panel of Figure 1, chromatography on Hypercarb column promotes additional higher-charge state precursor ions than C₁₈ which is instrumental in obtaining a high quality ETD MS/MS. Overall intensity of glycopeptide precursor ions with Hypercarb column are the same or higher than with C₁₈ column.

