

# Measuring Yield Stress to Correlate Slump of Concrete and Cement Paste

Jan Philip Plog, Thermo Fisher Scientific, Material Characterization, Karlsruhe, Germany

## Key words

Rheology, Yield Stress, Concrete, Cement Paste, Slump Test

## Introduction

Concrete is the most widely used man-made material (measured by tonnage made) in the world. One of the most important characteristics of these materials is the so-called workability. The workability is the ability of a fresh (plastic) concrete mix to fill a given form/mold properly with the desired work (vibration) and without reducing the concrete's quality. Workability depends amongst others on water content, aggregate (shape and size distribution), cementitious content and age (level of hydration) and can be modified by adding chemical admixtures, like super plasticizer or raising water content. However, excessive water leads to increased bleeding (build-up of surface water) and/or segregation aggregates (when the cement and aggregates start to separate), with the resulting concrete having reduced quality [1].

Workability can be measured by the "concrete slump test", a simplistic measure of the plasticity of a fresh batch of concrete following ASTM C 143 [2] or EN 12350-2 [3] test standards. Slump is normally measured by filling an "Abrams cone" with a sample from a fresh batch of concrete. The cone is placed with the wide end down onto a level, non-absorptive surface. It is then filled in three layers of equal volume, with each layer being tamped with a steel rod to consolidate the layer. When the cone is carefully lifted off, the enclosed material slumps a certain amount, owing to gravity. However, it was shown in the past that determining the yield stress rheologically is a fast and easy method to correlate slump [4]. As soft foods like concrete are often difficult to work with when using conventional plate/plate or concentric cylinder geometries on rotational rheometers because of the possible wall slip and excessive sample disruption during loading into narrow gaps vane geometries are recommended here.

When a vane rotor is fully immersed in the sample, the yield stress itself can then be calculated according to Boger [5]:

$$\tau_0 = \frac{M_{\max}}{K} \quad [a]$$

With M being the maximum Torque and K the vane parameter that depends on the height (H) and the diameter (D)



Fig. 1: Thermo Scientific HAAKE Viscotester iQ

of the paddle according to:

$$K = \frac{\pi \cdot D^3}{2} \left[ \frac{H}{D} + \frac{1}{3} \right] \quad [b]$$

## Experimental Results and Discussion

As described earlier it is recommendable to rheologically test concrete with vane rotors. Figure 1 shows the new Thermo Scientific™ HAAKE™ Viscotester™ iQ with vane configuration.

For the tests, a standard Portland cement was mixed with water at typical concentrations and then fine gravel was added at three different concentrations w/w. Please find an overview of the sample compositions in Table 1.

Sample	Portland Cement/g	Water/g	Fine Gravel/g
1	275	75	0
2	275	75	125
3	275	75	75

Table 1: Composition of tested concrete samples.




Fig. 2: Shear Stress versus Time for the three different concrete formulations at 25 °C after 5 min. waiting period after mixing.

After the samples have been formulated by thorough mixing they have been tested after 5 min. waiting time. The test was conducted as follows.

After the vane rotor has been fully immersed into the sample a constant rotational speed  $\Omega = 0.05$  rpm was commanded. Then the shear stress is monitored as a function of measuring time. After an initial purely elastic response in the sample the structure collapses and the shear stress is decreasing again. The maximum value in shear stress then corresponds with the yield stress. Figure 2 shows the results for the three formulations 5 min. after mixing at 25 °C.

As can be seen in Figure 2 the yield stress for the pure cement paste is 375 Pa compared to 890 Pa with 75 g of fine gravel and 1630 Pa with 125 g of fine gravel. Those easily determined yield stress values can now be transferred into the slump values (mm) determined with the ASTM Abrams cone according to the semi-empirical relationship by Hu et. al. [6]:

$$s = 300 - 347 \frac{(\tau_0 - 212)}{\rho} \quad [c]$$

with  $s$  being the slump in mm,  $\tau_0$  being the yield stress and  $\rho$  being the density of the fluid.

## Conclusion

The vane method on the Thermo Scientific HAAKE Viscotester iQ is a quick, simple and accurate approach to measure the yield stress of cement paste and concrete. Those values can then be easily transferred into slump values via semi-empirical relationships.

## Literature

- [1] See Wikipedia, Concrete, <http://en.wikipedia/wiki/Concrete> (as of Jan. 13, 2014, 16:50 CET)
- [2] ASTM C 143 - Standard Test Method for Slump of Hydraulic-Cement Concrete
- [3] EN 12350-2 - Testing fresh concrete. Slump test
- [4] N. Roussel, Correlation between Yield Stress and Slump: Comparison between Numerical Simulations and Concrete Rheometers Results, Materials and Structures, May 2006, Volume 39, Issue 4, pp 501-509
- [5] Dzuy NQ, Boger DV. 1985. Direct yield stress measurement with the vane method. J Rheol 29:335-47
- [6] Chong Hu, François De Larrard, Odd E. GjØrv, Rheological testing and modelling of fresh high performance concrete, Materials and Structures, January/February 1995, Volume 28, Issue 1, pp 1-7

[thermoscientific.com/mc](http://thermoscientific.com/mc)

© 2014/01 Thermo Fisher Scientific Inc. Copyrights in and to all photographs of instruments are owned by Thermo Fisher Scientific. This document is for informational purposes only. Specifications, terms and pricing are subject to change. Not all products are available in every country. Please consult your local sales representative for details.

### Material Characterization

**Benelux**  
Tel. +31 (0) 76 579 55 55  
info.mc.nl@thermofisher.com

**China**  
Tel. +86 (221) 68 65 45 88  
info.mc.china@thermofisher.com

**France**  
Tel. +33 (0) 1 60 92 48 00  
info.mc.fr@thermofisher.com

**India**  
Tel. +91 (20) 6626 7000  
info.mc.in@thermofisher.com

**Japan**  
Tel. +81 (45) 453-9167  
info.mc.jp@thermofisher.com

**United Kingdom**  
Tel. +44 (0) 1606 548 100  
info.mc.uk@thermofisher.com

**USA**  
Tel. +1 603 436 9444  
info.mc.us@thermofisher.com

**International/Germany**  
Dieselstr.4  
76227 Karlsruhe  
Tel. +49 (0) 721 4 09 44 44  
info.mc.de@thermofisher.com

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific