



The following features are present in the vector based on the nucleotide sequence.

Cytomegalovirus (CMV) promoter.....	9-636
Red Firefly luciferase gene	647-2293
BGH poly(A) signal.....	2305-2430
SV40 origin/promoter.....	2431-2995
Puromycin resistant gene	2996-3595
SV40 poly(A) signal	3757-3790
beta-lactamase (Amp ^R) gene	3899-4759
pUC replication origin (pUC Ori).....	4938-5742
Transcriptional terminator (Ter)	5743-6350
Lac operator 1 (Lac O1)	6351-6371
Transcriptional pause site (TPS).....	6504-6575