



The following features are present in the vector based on the nucleotide sequence.

CMV promoter	18-618
T7 promoter	625-642
Tluc.....	701-1195
Bovine growth hormone (BGH) poly(A) signal	1231-1356
f1 origin	1546-1852
Simian virus 40 (SV40) origin/promoter.....	2100-2359
Kanamycin/neomycin resistance gene.....	2394-3188
Herpes simplex virus (HSV) thymidine kinase (TK) poly(A) signal	3424-3442
pUC origin.....	3773-4416