



The following features are present in the vector based on the nucleotide sequence.

CMV promoter .....	18-618
T7 promoter .....	625-642
Tluc-Bcl-xL .....	665-1891
Bovine growth hormone (BGH) poly(A) signal .....	1927-2052
f1 origin .....	2242-2548
Simian virus 40 (SV40) origin/promoter .....	2796-3055
Kanamycin/neomycin resistance gene .....	3090-3884
Herpes simplex virus (HSV) thymidine kinase (TK) poly(A) signal .....	4120-4138
pUC origin .....	4469-5112