



Comments for pBlueBac4.5:

4940 nucleotides

- Polyhedrin promoter (P_{PH}): bases 7-95
- Multiple cloning site: bases 132-198
- SV40 polyadenylation sequence: 199-326
- Recombination sequences (ORF1629): bases 1343-543
- Ampicillin resistance gene: bases 1730-2587
- ColE1 origin: bases 2735-3408
- 5' lacZ fragment: bases 4638-3526
- lacZ sequence homologous to lacZ sequence
in Bac-N-Blue™ DNA: bases 4414-3526
- Early-to-late promoter (P_{ETL}): bases 4940-4639