

Comments for pBlueBac4.5/CAT:

5588 nucleotides



Polyhedrin promoter (P_{PH}): bases 7-95

CAT ORF: bases 222-878

SV40 polyadenylation sequence: 986-1113

Recombination sequences (ORF1629): bases 2130-1330

Ampicillin resistance gene: bases 2517-3374


ColE1 origin: bases 3522-4195

5' *lacZ* fragment: bases 5425-4313

lacZ sequence homologous to *lacZ* sequence

in Bac-N-Blue™ DNA: bases 5201-4313

Early-to-late promoter (P_{ETL}): bases 5727-5426



24-0278 A-160812

The sequence of pBlueBac4.5/CAT has been compiled from information in sequence databases, published sequences, and other sources. This vector has not been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Services Department.