

Comments for pBlueBac4/CAT:

5588 nucleotides



Polyhedrin promoter (P_{PH}): bases 7-95

ATG initiation codon: bases 168-170

CAT ORF: bases 217-873

Recombination sequences (ORF1629): bases 1991-1191


Ampicillin resistance gene: bases 2378-3235

ColE1 origin: bases 3383-4056

5' *lacZ* fragment: Bases 5286-4174

lacZ sequence homologous to *lacZ* sequence
in Bac-N-Blue™ DNA: bases 5062-4174

Early-to-late promoter (P_{ETL}): bases 5588-5287



vB-150721

The sequence of pBlueBac4/CAT has been compiled from information in sequence databases, published sequences, and other sources. This vector has not yet been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Services Department at 800-955-6288.