


**Comments for pCDM8:  
4356 base pairs**



pMB1 origin (pBR322): bases 1-587  
M13 origin: bases 588-1182  
SupF gene: bases 1183-1384  
Synthetic supF gene: bases 1183-1384  
CMV promoter: bases 1516-2170  
T7 promoter: bases 2171-2189  
Polylinker and stuffer fragment: bases 2198-2655  
Splice and polyadenylation signals: bases 2656-3352  
Polyoma origin and enhancers: bases 3353-4188  
SV40 origin (enhancerless): bases 4189-4352



vA-151005

The sequence of pCDM8 has been compiled from information in sequence databases, published sequences, and other sources. Portions of this vector have not yet been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Services Department at 800-955-6288.