


**Comments for pcDNA1.1/Amp:**  
4767 nucleotides



Col E1 origin: bases 1-587  
M13 origin: bases 588-1180  
Ampicillin gene: bases 1377-2235  
CMV promoter: bases 2304-2954  
T7 primer sequence: bases 2938-2957  
Polylinker: bases 2964-3063  
Splice and PolyA: bases 3060-3758  
Polyoma origin: bases 3764-4659  
SV40 origin: bases 4600-4767



B-180221

The sequence of pcDNA1.1/Amp has been compiled from information in sequence databases, published sequences, and other sources. This vector has not been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Service Department.