

Comments for pcDNA™3.1/Zeo (+)
5015 base pairs

CMV promoter: bases 209-863
 T7 promoter/priming site: bases 863-882
 Multiple cloning site: bases 895-1010
 BGH reverse priming site: bases 1022-1039
 BGH polyadenylation signal: bases 1021-1235
 f1 origin: bases 1298-1711
 SV40 promoter and origin: bases 1776-2101
 EM7 promoter: bases 2117-2183
 Zeocin™ resistance gene: bases 2184-2558
 SV40 polyadenylation signal: bases 2688-2817
 pUC origin: bases 3201-3874 (c)
 bla promoter: bases 4880-4978 (c)
 Ampicillin resistance (bla) gene: bases 4019-4879 (c)

(c) = complementary strand

