



CMV promoter: bases 232-819

T7 promoter/priming site: bases 863-883

*loxH* site: bases 901-934

BGH reverse priming site: bases 954-971

BGH polyadenylation sequence: bases 960-1184

f1 origin: bases 1230-1658

SV40 early promoter and origin: bases 1683-1992

EM-7 promoter: bases 2040-2095

Zeocin resistance gene: bases 2114-2488

SV40 late polyadenylation signal: bases 2618-2749

pUC origin: bases 3046-3949

*bla* promoter: bases 4810-4908 (C)

Ampicillin (*bla*) resistance gene: bases 3949-4809 (C)

(C) = complementary strand