



Comments for pcDNA™6.2/GFP-GW/p64^{TAG}:
7103 nucleotides

CMV promoter: bases 232-819
 T7 promoter/priming site: bases 863-882
attB1 site: bases 911-935
 p64^{TAG} ORF: bases 938-2257
attB2 site: bases 2259-2283
 Cycle-3 GFP: bases 2303-3016
 GFP reverse priming site: 2411-2432
 TK polyadenylation signal: bases 3031-3302
 f1 origin: 3338-3766
 SV40 early promoter and origin: 3793-4101
 EM7 promoter: bases 4156-4222
 Blasticidin resistance gene: bases 4223-4621
 SV40 early polyadenylation signal: bases 4779-4909
 PUC origin: bases 5292-5962 (c)
 Ampicillin (*bla*) resistance gene: bases 6107-6967 (c)
bla promoter: bases 6968-7066 (c)

(c) = complementary strand