



Comments for pcDNA6/TR[®]
6662 nucleotides

CMV promoter: bases 232-819

Rabbit β-globin intron II (IVS): bases 1028-1600

TetR gene: bases 1684-2340

SV40 early polyadenylation sequence: bases 2346-2477

f1 origin: bases 2897-3325

SV40 promoter and origin: bases 3335-3675

EM-7 promoter: bases 3715-3781

Blasticidin resistance gene: bases 3782-4180

SV40 early polyadenylation sequence: bases 4338-4468

pUC origin: bases 4851-5521

bla promoter: bases 6521-6625 (complementary strand)

Ampicillin (*bla*) resistance gene: bases 5666-6526 (complementary strand)