


**Comments for pcDNA II:  
3013 nucleotides**



Polylinker: bases 10-122  
Sp6 promoter: bases 136-152  
Ampicillin resistance: bases 1331-2191  
F1 origin: bases 2377-2832  
Lac Z gene: bases 2832-390  
T7 promoter: bases 2993-3012



vA-951003

The sequence of pcDNA II has been compiled from information in sequence databases, published sequences, and other sources. Portions of this vector have not yet been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Services Department at 800-955-6288.

**U.S. Headquarters**

Tel: 1-800-955-6288  
Fax: 1-619-597-6201

**European Headquarters**

Tel: +31 (0) 5945-15175  
Fax: +31 (0) 5945-15312