



Comments for pCEP4/CAT:
10891 nucleotides

CMV promoter: bases 1-588

Chloramphenicol acetyl transferase (CAT) gene: bases 675-1334

SV40 polyadenylation signal: bases 1390-1631

OriP: bases 2049-4024

EBNA-1 gene (complementary strand): bases 4325-6250

Ampicillin resistance gene: bases 6876-7736

pUC origin: bases 7745-8520

TK promoter: bases 8888-9050

Hygromycin resistance gene: bases 9114-10124

TK polyadenylation signal: bases 10136-10407