



**Comments for:  
(no insert)**

T7 promoter/priming site:  
 Sp6 promoter/priming site (c):  
 CMV promoter (c):  
*lacUV5* promoter (c):  
 pUC origin:  
 Ampicillin (*bla*) resistance gene (c):  
*bla* promoter (c):  
 f1 origin (c):  
 SV40 early polyadenylation signal (c):

**pCMV•SPORT1  
4380 nucleotides**

16-35  
 156-175  
 262-855  
 856-883  
 1123-1821  
 2146-3006  
 3001-3058  
 3294-3749  
 3876-4142

**pCMV•SPORT2  
4475 nucleotides**

11-30  
 151-170  
 252-845  
 900-927  
 1335-2004  
 2191-3051  
 3046-3103  
 3339-3794  
 3921-4187

**pCMV•SPORT4  
4387 nucleotides**

11-30  
 151-169  
 252-845  
 --  
 1243-1916  
 2103-2963  
 2959-3015  
 3251-3706  
 3833-4099

(c) = complementary strand