



**Comments for pCMV•SPORT6.1 (no insert)  
4177 nucleotides**

SV40 early polyadenylation signal (c): bases 126-392

M13 Forward (-20) priming site: bases 664-679

T7 promoter/priming site: bases 691-710

*attB2*: bases 711-735

*attB1*: bases 807-831

Sp6 promoter/priming site (c): bases 833-850

M13 Reverse priming site (c): bases 861-877

CMV promoter (c): bases 920-1510

pUC origin: bases 1895-2562

Ampicillin (*bla*) resistance gene (c): bases 2729-3589

*bla* promoter (c): bases 3590-3641

f1 origin (c): bases 3721-4176

(c) = complementary strand