



**Comments for pCR®2.1/GAL4-VP16**  
**5842 nucleotides**

*lac* promoter: bases 143-172  
 SV40 early promoter: bases 360-660  
 GAL4 DNA binding domain: bases 676-1116  
 VP16 activation domain: bases 1135-1368  
 SV40 early polyadenylation sequence: bases 1898-2141  
 LacZ $\alpha$  fragment: bases 2321-2481  
 f1 origin: bases 2482-2896  
 Kanamycin promoter: bases 3092-3229  
 Kanamycin resistance ORF: bases 3230-4024  
*bla* promoter: bases 4025-4047  
 Ampicillin resistance ORF: bases 4042-4092  
 pUC origin: bases 5047-5720