



Comments for pCR[®]2.1/VP16-CP
5504 nucleotides

lac promoter: bases 143-172

SV40 early polyadenylation sequence: bases 389-519 (C)

Polyoma viral coat protein, VP3 (amino acids 1-204): bases 631-1242 (C)

VP16 activation domain: bases 1243-1509 (C)

SV40 early promoter: bases 1537-1858 (C)

LacZ α fragment: bases 1983-2143

f1 origin: bases 2144-2558

Kanamycin promoter: bases 2754-2891

Kanamycin resistance ORF: bases 2892-3686

bla promoter: bases 3687-3710

Ampicillin resistance ORF: bases 3704-4564

pUC origin: bases 4709-5384

(C)= complementary strand