



CMV promoter: bases 1-596

Putative transcriptional start: bases 620-625

T7 promoter/priming site: bases 638-657

CAT ORF: bases 738-1486

BGH reverse priming site: bases 1562-1579

BGH polyadenylation site: bases 1561-1775

pUC origin: bases 1906-2489

SV40 promoter and origin: bases 3943-4281 (complement)

Neomycin/kanamycin resistance gene (ORF): bases 3120-3908 (complement)

Thymidine kinase polyadenylation site: bases 2675-2945 (complement)

Ampicillin resistance gene (ORF): bases 4360-5220 (complement)

f1 origin: bases 5351-5807