

Comments for pDEST™15
7013 base pairs

T7 promoter: bases 25-44

Ribosome binding site (RBS): bases 90-96

Initiation ATG: bases 105-107

GST tag: bases 108-776

attR1: bases 792-916

Chloramphenicol resistance gene (*Cm^R*): bases 1025-1684

ccdB gene: bases 2026-2331

attR2: bases 2372-2496

T7 transcription termination region: bases 2518-2646

bla promoter: bases 3134-3232

Ampicillin resistance (*bla*) gene: bases 3233-4093

pBR322 origin: bases 4238-4911

rop ORF: bases 5282-5473 (c)

(c) = complementary strand

