



**Comments for pEF6/Myc-His/lacZ
8980 nucleotides**

- EF-1 α promoter: bases 388-1547
 - T7 promoter/priming site: bases 1548-1567
 - LacZ ORF: bases 1648-4704
 - myc epitope: bases 4729-4758
 - Polyhistidine tag: bases 4774-4791
 - BGH reverse priming site: bases 4814-4831
 - BGH polyadenylation signal: bases 4820-5044
 - f1 origin: bases 5090-5518
 - SV40 promoter and origin: bases 5545-5853
 - EM-7 promoter: bases 5908-5974
 - Blasticidin resistance gene: bases 5975-6373
 - SV40 polyadenylation signal: bases 6531-6661
 - pUC origin: bases 7044-7716 (C)
 - b/a promoter: bases 8757-8841 (C)
 - Ampicillin resistance gene: bases 7862-8722 (C)
- (C) = complementary strand