



**Comments for pEF/Bsd:
4304 nucleotides**

- bla* promoter: bases 102-200
- Ampicillin (*bla*) resistance gene: bases 201-1061
- pUC origin: bases 1206-1879
- 5' polylinker: bases 1884-1928
- Human EF-1 α promoter: bases 2020-3198
- EM7 promoter: bases 3222-3288
- Blasticidin (*bsd*) resistance gene: bases 3289-3687
- SV40 early polyadenylation sequence: bases 3729-3859
- 3' polylinker: bases 3859-3904