



EF-1 α promoter: bases 6-1192

EF Forward priming site: bases 1133-1153

GFP ORF: bases 1220-1936

myc epitope: bases 1946-1975

BGH Reverse priming site: 2011-2028

BGH polyadenylation sequence: bases 2010-2224

f1 origin: bases 2287-2700

SV40 promoter/origin: bases 2765-3072

Neomycin resistance gene (ORF): bases 3080-3874

SV40 polyadenylation sequence: bases 3891-4129

pUC origin: bases 4312-4985 (Complementary strand)

Ampicillin resistance gene (ORF): 5130-5990 (Complementary strand)