



EF-1 α promoter: bases 6-1192

EF Forward priming site: bases 1133-1153

Mitochondrial targeting sequence: bases 1196-1282

GFP ORF: bases 1289-2005

myc epitope: bases 2015-2044

BGH Reverse priming site: 2080-2097

BGH polyadenylation sequence: bases 2079-2293

f1 origin: bases 2356-2769

SV40 promoter/origin: bases 2834-3141

Neomycin resistance gene (ORF): bases 3149-3943

SV40 polyadenylation sequence: bases 3960-4198

pUC origin: bases 4381-5054 (Complementary strand)

Ampicillin resistance gene (ORF): 5199-6059 (Complementary strand)