



**Comments for pEXP4-DEST
4415 nucleotides**

- T7 promoter: bases 20-36
- T7 promoter/priming site: bases 20-39
- attR1 site: bases 98-222
- Chloramphenicol resistance gene (Cm^R): bases 331-990
- ccdB gene: bases 1332-1637
- attR2 site: bases 1678-1802
- Lumio™ tag: bases 1831-1848
- Polyhistidine (6xHis) region: bases 1864-1881
- T7 reverse priming site: bases 1935-1954
- T7 transcription termination region: bases 1940-2024
- Zeocin™ resistance gene: bases 2168-2542
- Ampicillin resistance gene: bases 2563-3423
- pUC origin: bases 3568-4241