



Comments for pFRT/lacZeo2
6482 nucleotides

SV40 Δ early promoter and origin: bases 36-248
 ATG initiation codon: bases 305-307
 FRT site: bases 310-357
LacZ-Zeocin™ fusion gene
LacZ ORF (no ATG): bases 374-3418
 Zeocin™ resistance gene (no ATG): bases 3506-3877
 SV40 early polyadenylation signal: bases 4007-4137
bla promoter: bases 4684-4782
Ampicillin (bla) resistance gene: bases 4783-5643
 pUC origin: bases 5788-6461