



**Comments for pIND/Hygro
5193 nucleotides**

Ecdysone/glucocorticoid Response Elements (5x E/GREs): bases 12-185
 Minimal Heat Shock Promoter and 5' UTR: bases 192-486
 Putative Transcriptional Start: base 245
 Multiple Cloning Site: bases 487-612
 BGH Reverse priming site: bases 614-631
 BGH Polyadenylation sequence: bases 613-827
 f1 origin: bases 890-1303
 SV40 promoter/origin: bases 1368-1693
 Hygromycin-B-phosphotransferase ORF (Hyg^R): bases 1711-2734
 SV40 Polyadenylation sequence: bases 2747-3119
 pUC origin: bases 3379-4052 (Complementary strand)
 β-lactamase ORF (Amp^R): bases 4197-5057 (Complementary strand)