



Comments for pIND/Hygro/lacZ
8244 nucleotides

Ecdysone/glucocorticoid

Response Elements (5x E/GREs): bases 12-185

Minimal Heat Shock Promoter and 5' UTR: bases 192-486

Putative Transcriptional Start: base 245

LacZ ORF: bases 564-3636

BGH Reverse priming site: bases 3665-3682

BGH Polyadenylation sequence: bases 3664-3878

f1 origin: bases 3941-4354

SV40 promoter/origin: bases 4419-4744

Hygromycin-B-phosphotransferase

ORF (Hyg^R): bases 4762-5785

SV40 Polyadenylation sequence: bases 5798-6170

pUC origin: bases 6430-7103 (C)

β-lactamase ORF (Amp^R): bases 7248-8108 (C)

(C) = complementary strand