



**Comments for pIZ/V5-His
2876 nucleotides**

OpIE2 promoter: bases 4-552
 Multiple cloning site: bases 561-656
 V5 epitope: bases 663-704
 6xHis tag: bases 714-731
 OpIE2 Reverse priming site: bases 741-766
 OpIE2 polyadenylation sequence: bases 749-878
 pUC origin: bases 947-1620
 OpIE2 promoter: bases 1665-2213
 EM7 promoter: bases 2231-2308
 Zeocin™ resistance marker (ORF): bases 2309-2683
 SV40 early polyadenylation sequence: bases 2747-2876