



Comments for pPIC6 A:
3382 nucleotides

5' AOX1 promoter region: bases 1-942
 5' AOX1 priming site: bases 855-875
 Multiple cloning site: bases 932-1011
c-myc epitope: bases 1012-1041
 Polyhistidine (6xHis) tag: bases 1057-1074
 3' AOX1 priming site: bases 1160-1180
 AOX1 transcription termination region: bases 1078-1419
 TEF1 promoter: bases 1420-1828
 EM7 promoter: bases 1833-1899
 Blasticidin resistance gene: bases 1900-2298
 CYC1 transcription termination region: bases 2327-2644
 pUC origin: bases 2655-3328 (complementary strand)

* The restriction site between *Not* I and the *c-myc* epitope is different in each version of pPIC6:

Apa I in pPIC6 A
Xba I in pPIC6 B
SnaB I in pPIC6 C