



Comments for pPICZ-E/Uni-*lacZ*
8550 nucleotides

5' AOX1 promoter region: bases 1-942
 5' AOX1 priming site: bases 855-875
loxP site: bases 949-982
lacZ ORF: bases 1003-4062
 BGH polyadenylation region: bases 4222-4430
 T7 transcriptional termination region: bases 4445-4573
 Kanamycin resistance gene: bases 4752-5546 (C)
 kan promoter region: bases 5547-5684 (C)
 R6K γ origin: bases 6018-6290
loxP site: bases 6319-6352
 AOX1 transcription termination region: bases 6353-6639
TEF1 promoter: bases 6640-7048
 EM7 promoter: bases 7053-7101
Zeocin resistance gene: bases 7120-7494
 CYC1 transcription termination region: bases 7495-7812
 pUC origin: bases 7822-8496
 (C) =complementary strand