



**Comments for pREP7:
9283 nucleotides**

- RSV LTR promoter: bases 58-628
- Multiple cloning site: bases 632-691
- SV40 polyadenylation signal: bases 684-925
- TK polyadenylation signal: bases 1332-1603
- Hygromycin resistance gene (complementary strand): bases 1615-2625
- TK promoter: bases 2689-2851
- pUC origin: bases 3219-3994
- Ampicillin resistance gene (complementary strand): bases 4003-4863
- EBNA-1 gene: bases 5489-6739
- OriP: bases 7043-9018