


**Comments for pSE420:
4617 nucleotides**



Trc promoter: bases 191-221
 Lac operator: bases 228-248
 rrnB antitermination sequences: bases 264-333
 T7 gene 10 translational enhancer: bases 346-354
 Ribosome binding site: bases 370-374
 Mini-cistron: bases 383-409
 Superlinker MCS: bases 411-763
 Ampicillin resistance ORF: bases 1283-2143
 LacI^q ORF: bases 3615-4574



The sequence of pSE420 has been compiled from information in sequence databases, published sequences, and other sources. Portions of this vector have not yet been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Services Department at 800-955-6288.

**For Technical Service Call
1-800-955-6288**