

Comments for pSecTag A:

5166 nucleotides

CMV promoter: bases 209-863

T7 promoter/priming site: bases 863-882

Ig κ -chain leader sequence: bases 905-967

Multiple cloning site: bases 968-1088

myc epitope: bases 1089-1121

Polyhistidine tag: bases 1134-1151

pcDNA3.1/BGH reverse priming site: bases 1174-1191

BGH poly A: bases 1173-1387

f1 origin of replication: bases 1298-1711


SV40 promoter and origin: bases 1931-2252

ZeocinTM resistance gene: bases 2335-2709

SV40 poly A: bases 2839-2968

ColE1 origin: bases 3352-4025

Ampicillin resistance gene: bases 4170-5030



A-160422

The sequences of pSecTag A, B, and C have been compiled from information in sequence databases, published sequences, and other sources. These vectors have not been completely sequenced. If you suspect an error in the sequences, please contact Invitrogen's Technical Services Department.