

Comments for pSinRep5:

9951 nucleotides

SP6 promoter: bases 9933-9951

Nonstructural genes: bases 60-7598

Subgenomic promoter (P_{SG}): bases 7580-7603

Forward Sindbis sequencing priming site: bases 7604-7623

Transcriptional start: base 7598


Multiple cloning site: bases 7647-7689

PolyA tail: bases 7997-8033

3' Restriction sites: bases 8047-8067

Ampicillin resistance gene: bases 8227-9085

ColE1 origin: bases 9232-9861



C-180402

The sequence of pSinRep5 has been compiled from information in sequence databases, published sequences, and other sources. This vector has not yet been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Services Department.