



trc promoter: bases 191-221
lac operator: bases 228-248
rrnB anti-termination sequences: bases 264-333
 T7 gene 10 translational enhancer: bases 346-354
 Ribosome binding site: bases 370-374
 Mini-cistron: bases 383-409
 Polyhistidine and enterokinase cleavage site: bases 425-504
 Xpress™ epitope: bases 482-505
 Multiple cloning site: bases 515-554
rrnB transcriptional termination sequence: bases 637-794
 Ampicillin resistance ORF: bases 1074-1934
 pBR322 origin: bases 2079-2752
lac I^q ORF: bases 3281-4365