


Comments for pVL1392:
9632 nucleotides

Recombination sequence (ORF603+): bases 1-3997
Polyhedrin promoter: bases 3998-4092
Polyhedrin gene: bases 4093-4738
Multiple cloning site: bases 4128-4179
Recombination sequence (ORF1629+): bases 4738-7002
ColE1 origin: bases 8029-7356
Ampicillin resistance gene: bases 8965-8177



A-150622

The sequence of pVL1392 has been compiled from information in sequence databases, published sequences, and other sources. Portions of this vector have not yet been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Services Department at 800-955-6288.