



Comments for pZeoSV2
3515 nucleotides

- SV40 promoter: bases 10-356
- T7 promoter/primer: bases 401-420
- Multiple Cloning Site: bases 433-548
- BGH reverse priming site: bases 560-577
- BGH polyA: bases 560-773
- pUC origin: bases 814-1487
- f1 origin: bases 1780-2290
- CMV promoter: bases 2304-2930
- EM-7 promoter: bases 2931-2997
- Zeocin resistance gene: 2998-3372
- SV40 polyA: bases 3380-3515